

TA IQ


COMPANY:		Shipping address:	
Customer number:	Date:	Address:	
		Postcode:	Town/City:
Cost Estimate:	Order:	Country:	Phone:
Contact person:		E-Mail:	
Reference:			

Specifications

	FWD	MWD	RWD
Lenght with leg support: (mm)	985-1100	985-1100	940-1060
Overall width: (mm)	630 - 735 for LC		
Seat width min.-max.: (mm)	250-550		
Seat depht min.-max.: (mm)	250-600		
Seat height w/o cushion: (mm)	380		
Backrest height min.-max.: (mm)	380-790		
Seat tilt: (°)	45		
Backrest angle min.-max.: (°)	85-135		
Lift: (mm)	300		
Footplate to seat distance min.-max.: (mm)	110-530		
Leg to seat surface angle min.-max.: (°)	90-180		
Arm support to seat distance min.-max.: (mm)	370-475		

	FWD	MWD	RWD
Total mass: (kg)	150	155 -165	150
Maximum user weight: (kg)	140		
Dynamic stability uphill: (°)	10		
Obstacle climbing: (mm)	100	80	80
Ground clearance: (mm)	700		
Max. speed forward: (km/h)	10 / 12 / 15		
Minimum turning radius: (mm)	650	450	650

User details



All seat dimensions are always measured without seat cushion!

01 Seat depth:
02 Seat width:
03 Chest width:
04 Chest depth:
05 Elbow height:
06 Arm pit height:
07 Back height:
08 Head height:
09 Forearm length:
10 Lower leg length:
(Measured from seatplate to footplate)
Height:
Weight (kg):

Specify dimensions in mm
Refer to the sketches!

The grey colored rows are indicative and show a standard configuration

1 Chassis / Motor

	CODE	DESCRIPTION		Speed	Motor Power	Range (*)	Max. User weight	PRICE
☐	97050	TA IQ Midwheel electric wheelchair Lift 300 mm & tilt 45° included		12 km/h	2 x 350 W	40 km	140 kg	€ 3.371,00
☐	97055	TA IQ Midwheel electric wheelchair Lift 300 mm & tilt 45° included		10 km/h	2 x 350 W	40 km	140 kg	€ 3.371,00
☐	97056	TA IQ Midwheel LC electric wheelchair including drive wheels 4.00-6, and 4" black pneumatic caster Lift 300 mm & tilt 45° included		15 km/h	2 x 600 W	20-25 km ¹⁾ 30-35 km ²⁾	140 kg	€ 5.080,00
☐	97057	TA IQ Midwheel RS electric wheelchair Lift 300 mm & tilt 45° included		15 km/h	2 x 600 W	25-30 km ¹⁾ 35-40 km ²⁾	140 kg	€ 4.280,00
☐	97060	TA IQ Rearwheel electric wheelchair Lift 300 mm & tilt 45° included		12 km/h	2 x 350 W	40 km	140 kg	€ 3.370,00
☐	97065	TA IQ Rearwheel electric wheelchair Lift 300 mm & tilt 45° included		10 km/h	2 x 350 W	40 km	140 kg	€ 3.370,00
☐	97070	TA IQ Frontwheel electric wheelchair Lift 300 mm & tilt 45° included		12 km/h	2 x 350 W	40 km	140 kg	€ 3.370,00
☐	97075	TA IQ Frontwheel electric wheelchair Lift 300 mm & tilt 45° included		10 km/h	2 x 350 W	40 km	140 kg	€ 3.370,00

(*) According TA requirements



97050/97055/97057



97060/97065



97070/97075



97056

¹⁾ with 80 Ah Battery

²⁾ with 100 Ah Battery

Options for chassis / Motor

	CODE	DESCRIPTION	PRICE
☐	95380	Gyro module (standard incl. for FWD) includes item: 95286, 95265 and 95507	€ 288,00
☐	97349	Antitip for RWD and FWD, set	€ 230,00
☐	97266-1	Dahl car docking with lock plate V2	€ 350,00
☐	97267-1	Dahl car attachment w/o lock plate V2	€ 250,00



97266-1

2 Battery and chargers

	CODE	DESCRIPTION	PRICE
☐	97432	Batteries 80Ah GEL (2 pcs)	pcs. price € 200,00
☐	97322	Batteries 100Ah AGM(2 pcs)	pcs. price € 270,00
☐	97190	Charger 8A GEL	€ 190,00
☐	97284	Charger 10A Medico GEL	€ 240,00
☐	97374	Charger 24V-12A AGM	€ 330,00
☐	M1060475	Power cable Type J(for 8A charger)	Switzerland € 30,00
☐	M1060477	Power cable Type G (for 8A charger)	England, Ireland, Cyprus, Qatar, United Arab Emirates (U.A.E.), € 30,00
☐	M1060478	Power cable Type I (for 8A charger)	New Zealand, Australia, Argentina € 30,00

3 Wheels



Drive wheels 3.00-8
(not compatible with LC Version (97056))

	Pneumatic tyres		Puncture-proof tyres	
	Tire grey	Tire black	Tire grey	Tire black
Rims alu grey	☐ 97142 € 56,00	☐ 97263 € 56,00	☐ 97226 € 79,00	☐ 97264 € 79,00
Rims alu black	☐ 97244 € 56,00	☐ 97194 € 56,00	☐ 97265 € 79,00	☐ 97196 € 79,00

Castor wheels 200 x 50
(not compatible with LC Version (97056))

	Pneumatic tyres		Puncture-proof tyres	
	Tire grey	Tire black	Tire grey	Tire black
Rims alu grey	☐ 97292 € 33,00	☐ 97260 € 33,00	☐ 97227 € 45,00	☐ 97261 € 45,00
Rims alu black	☐ 97245 € 33,00	☐ 97259 € 33,00	☐ 97262 € 45,00	☐ 97195 € 45,00

Drive wheels 12½ x 2 1/4
(not compatible with LC Version (97056))

	Pneumatic tyres		Puncture-proof tyres	
	Tire grey	Tire black	Tire grey	Tire black
Rims alu black	☐ 97270 € 51,00	☐ 97131 € 51,00	☐ 97348 € 61,00	☐ 97143 € 61,00

Castor wheels 4"

	Pneumatic tyres		Puncture-proof tyres	
	Tire grey	Tire black	Tire grey	Tire black
Rims alu black	☐ 97269 € 51,00	☐ 97314 € 51,00	☐ 97313 € 61,00	☐ 97300 € 61,00

Options for chassis / Motor

CODE	DESCRIPTION	PRICE
97378	Standard fork, tire acts as protection (pair)	€ 0,00
97379		
☐ 97282	Guards castorwheels IQ small, set	Not in combination with wheels 4" € 36,00



97378/97379



97282

4 Seat

Lift / tilt

CODE	DESCRIPTION	PRICE
97104	Actuator lift 300 mm	€ 0,00
97099	Seat tilt 45°	€ 0,00
Electric seat glide 110mm		
		seat width
		380mm 20710-1
		400mm 20710-2
		450mm 20710-3
		500mm 20710-4
		550mm 20710-5
		€ 1.590,00



97104



97099

Seat width

CODE	DESCRIPTION	PRICE
☐	Seat width 250mm adjustable to 250 - 300 or 350 mm	€ 160,00
☐	Seat width 380mm including depth extender (80233)	€ 180,00
☐	Seat width 400mm including depth extender (80233)	€ 180,00
☐	Seat width 450mm including depth extender (80216)	€ 183,00
☐	Seat width 500mm including depth extender (80215) reinforced with middle plate (80502)	€ 220,00
☐	Seat width 550mm including depth extender (80215) reinforced with middle plate (80502)	€ 220,00

Seat depth

Seat depth is stepless adjustable from 250 - 600mm:

Specify required seat depth, mandatory: _____

Seat system

CODE	DESCRIPTION	Seat width					PRICE
☐	80519 TA Strap seat	380mm	400mm	450mm	500mm	550mm	€ 390,00
		☐	☐	☐	☐	☐	



80519

Seat cushion



Netti Sit

Good pressure distribution and stabilizing properties.
The foam properties are adapted to the size of the cushion – and therefore the user's weight.
No adaptation and adjustment necessary.
Temperature stabilizing cover available to secure a good micro climate and prevent pressure ulcers.

CODE	DESCRIPTION	WIDTH x DEPTH (mm)	PRICE
☐	N187306 Netti Sit 3D	380x400	€ 227,00
☐	N187307 Netti Sit 3D	380x450	€ 227,00
☐	N187308 Netti Sit 3D	380x500	€ 227,00
☐	N187309 Netti Sit 3D	400x400	€ 227,00
☐	N187310 Netti Sit 3D	400x450	€ 227,00
☐	N187311 Netti Sit 3D	400x500	€ 227,00
☐	N187333 Netti Sit 3D	400x550	€ 304,00
☐	N187312 Netti Sit 3D	430x400	€ 227,00
☐	N187313 Netti Sit 3D	430x450	€ 227,00
☐	N187314 Netti Sit 3D	430x500	€ 227,00
☐	N187315 Netti Sit 3D	450x400	€ 227,00
☐	N187316 Netti Sit 3D	450x450	€ 227,00
☐	N187317 Netti Sit 3D	450x500	€ 227,00
☐	N187318 Netti Sit 3D	450x550	€ 304,00
☐	N187319 Netti Sit 3D	500x450	€ 227,00
☐	N187320 Netti Sit 3D	500x500	€ 227,00
☐	N187321 Netti Sit 3D	500x550	€ 304,00
☐	N187322 Netti Sit 3D	550x500	€ 304,00
☐	N187323 Netti Sit 3D	550x550	€ 304,00
☐	N188306 Netti Sit Easy Care	380x400	€ 227,00
☐	N188307 Netti Sit Easy Care	380x450	€ 227,00
☐	N188308 Netti Sit Easy Care	380x500	€ 227,00
☐	N188309 Netti Sit Easy Care	400x400	€ 227,00
☐	N188310 Netti Sit Easy Care	400x450	€ 227,00
☐	N188311 Netti Sit Easy Care	400x500	€ 227,00
☐	N188333 Netti Sit Easy Care	400x550	€ 304,00
☐	N188312 Netti Sit Easy Care	430x400	€ 227,00
☐	N188313 Netti Sit Easy Care	430x450	€ 227,00
☐	N188314 Netti Sit Easy Care	430x500	€ 227,00
☐	N188315 Netti Sit Easy Care	450x400	€ 227,00
☐	N188316 Netti Sit Easy Care	450x450	€ 227,00
☐	N188317 Netti Sit Easy Care	450x500	€ 227,00
☐	N188318 Netti Sit Easy Care	450x550	€ 304,00
☐	N188319 Netti Sit Easy Care	500x450	€ 227,00
☐	N188320 Netti Sit Easy Care	500x500	€ 227,00
☐	N188321 Netti Sit Easy Care	500x550	€ 304,00
☐	N188322 Netti Sit Easy Care	550x500	€ 304,00
☐	N188323 Netti Sit Easy Care	550x550	€ 304,00



Varilite Evolution PSV

For wheelchair users with a high degree of tissue breakdown and symmetric positioning needs.

VARILITE's PSV takes guesswork out of cushion adjustment and makes it automatic. Evolution PSV uses VARILITE's revolutionary new Pressure Setting Valve (patent pending). The PSV works by sensing the amount of internal cushion pressure. As air is released, internal pressure decreases and the user is immersed in the cushion for effective pressure distribution.

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)	PRICE
☐	93000	Varilite Evolution PSV	300x350	€ 250,00
☐	93001	Varilite Evolution PSV	350x350	€ 250,00
☐	93002	Varilite Evolution PSV	350x400	€ 250,00
☐	93420	Varilite Evolution PSV	350x450	€ 250,00
☐	93421	Varilite Evolution PSV	350x500	€ 250,00
☐	93422	Varilite Evolution PSV	380x380	€ 250,00
☐	93423	Varilite Evolution PSV	380x430	€ 250,00
☐	93424	Varilite Evolution PSV	400x350	€ 250,00
☐	93003	Varilite Evolution PSV	400x400	€ 250,00
☐	93004	Varilite Evolution PSV	400x450	€ 250,00
☐	93005	Varilite Evolution PSV	400x500	€ 250,00
☐	93006	Varilite Evolution PSV	430x430	€ 250,00
☐	93007	Varilite Evolution PSV	450x400	€ 250,00
☐	93008	Varilite Evolution PSV	450x450	€ 250,00
☐	93009	Varilite Evolution PSV	450x500	€ 250,00
☐	93435	Varilite Evolution PSV	500x400	€ 250,00
☐	93268	Varilite Evolution PSV	500x450	€ 250,00
☐	93010	Varilite Evolution PSV	500x500	€ 250,00
☐	93277	Varilite Evolution PSV	550x450	€ 290,00
☐	93278	Varilite Evolution PSV	550x500	€ 290,00
☐	93279	Varilite Evolution PSV	600x450	€ 290,00
☐	93280	Varilite Evolution PSV	600x500	€ 290,00



Varilite Meridian TM

Therapists recommend the Meridian for individuals with motor and neurological dysfunction due to spinal-cord injury, stroke, multiple sclerosis, cerebral palsy, and traumatic brain injury. Therapists value the integral pre-ischial ridge and wedge effect created by the Meridian's dual chambers. User and caregivers value its low weight, comfort, and ease of use. No pumps or accessories. No kneading or manipulation.

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)	PRICE
☐	93362	Varilite Meridian TM	300x350	€ 250,00
☐	93363	Varilite Meridian TM	350x350	€ 250,00
☐	93364	Varilite Meridian TM	350x400	€ 250,00
☐	93411	Varilite Meridian TM	350x450	€ 250,00
☐	93412	Varilite Meridian TM	350x500	€ 250,00
☐	93414	Varilite Meridian TM	380x380	€ 250,00
☐	93415	Varilite Meridian TM	380x430	€ 250,00
☐	93413	Varilite Meridian TM	400x350	€ 250,00
☐	93365	Varilite Meridian TM	400x400	€ 250,00
☐	93366	Varilite Meridian TM	400x450	€ 250,00
☐	93367	Varilite Meridian TM	400x500	€ 250,00
☐	93368	Varilite Meridian TM	430x430	€ 250,00
☐	93369	Varilite Meridian TM	450x400	€ 250,00
☐	93370	Varilite Meridian TM	450x450	€ 250,00
☐	93371	Varilite Meridian TM	450x500	€ 250,00
☐	93571	Varilite Meridian TM	450x550	€ 250,00
☐	93410	Varilite Meridian TM	500x400	€ 250,00
☐	93372	Varilite Meridian TM	500x450	€ 250,00
☐	93373	Varilite Meridian TM	500x500	€ 250,00
☐	93416	Varilite Meridian TM	550x450	€ 290,00
☐	93417	Varilite Meridian TM	550x500	€ 290,00
☐	93418	Varilite Meridian TM	600x450	€ 290,00
☐	93419	Varilite Meridian TM	600x500	€ 290,00



Varilite Pro-Form NX™
Dual chamber

For wheelchair users with challenging, asymmetric skin protection, positioning and postural support requirements. Modifiable.

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)	PRICE
☐	93296	Varilite Pro-Form NX™, dual chamber	300x350	€ 290,00
☐	93297	Varilite Pro-Form NX™, dual chamber	350x350	€ 290,00
☐	93298	Varilite Pro-Form NX™, dual chamber	350x400	€ 290,00
☐	93113	Varilite Pro-Form NX™, dual chamber	380x380	€ 290,00
☐	93114	Varilite Pro-Form NX™, dual chamber	380x430	€ 290,00
☐	93299	Varilite Pro-Form NX™, dual chamber	400x400	€ 290,00
☐	93300	Varilite Pro-Form NX™, dual chamber	400x450	€ 290,00
☐	93301	Varilite Pro-Form NX™, dual chamber	400x500	€ 290,00
☐	93115	Varilite Pro-Form NX™, dual chamber	430x430	€ 290,00
☐	93302	Varilite Pro-Form NX™, dual chamber	450x400	€ 290,00
☐	93303	Varilite Pro-Form NX™, dual chamber	450x450	€ 290,00
☐	93304	Varilite Pro-Form NX™, dual chamber	450x500	€ 290,00
☐	93116	Varilite Pro-Form NX™, dual chamber	500x400	€ 290,00
☐	93117	Varilite Pro-Form NX™, dual chamber	500x450	€ 290,00
☐	93305	Varilite Pro-Form NX™, dual chamber	500x500	€ 290,00



Varilite Reflex™

Non-adjustable skin protection cushion featuring VARILITE's unique Air-Release Device™.

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)	PRICE
☐	93491	Varilite Reflex™	300x300	€ 130,00
☐	93492	Varilite Reflex™	300x350	€ 130,00
☐	93493	Varilite Reflex™	350x350	€ 130,00
☐	93494	Varilite Reflex™	350x400	€ 130,00
☐	93495	Varilite Reflex™	350x450	€ 130,00
☐	93496	Varilite Reflex™	350x500	€ 130,00
☐	93497	Varilite Reflex™	380x380	€ 130,00
☐	93498	Varilite Reflex™	380x430	€ 130,00
☐	93499	Varilite Reflex™	400x350	€ 130,00
☐	93500	Varilite Reflex™	400x400	€ 130,00
☐	93501	Varilite Reflex™	400x450	€ 130,00
☐	93502	Varilite Reflex™	400x500	€ 130,00
☐	93503	Varilite Reflex™	430x430	€ 130,00
☐	93504	Varilite Reflex™	450x400	€ 130,00
☐	93505	Varilite Reflex™	450x450	€ 130,00
☐	93506	Varilite Reflex™	450x500	€ 130,00
☐	93507	Varilite Reflex™	500x400	€ 130,00
☐	93508	Varilite Reflex™	500x450	€ 130,00
☐	93509	Varilite Reflex™	500x500	€ 130,00



Varilite Stratus

The Varilite Stratus Cushion is a great choice for people who have a moderate to low risk of tissue breakdown and moderate positioning needs. The cushion has a unique self-inflating valve, which allows air in and out of the cushion so it can contour to the user's shape.

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)	PRICE
☐	93047	Varilite Stratus	300x300	€ 140,00
☐	93048	Varilite Stratus	300x350	€ 140,00
☐	93049	Varilite Stratus	350x350	€ 140,00
☐	93050	Varilite Stratus	350x400	€ 140,00
☐	93110	Varilite Stratus	380x380	€ 140,00
☐	93111	Varilite Stratus	380x430	€ 140,00
☐	93051	Varilite Stratus	400x400	€ 140,00
☐	93052	Varilite Stratus	400x450	€ 140,00
☐	93053	Varilite Stratus	430x430	€ 140,00
☐	93054	Varilite Stratus	450x400	€ 140,00
☐	93055	Varilite Stratus	450x450	€ 140,00
☐	93112	Varilite Stratus	450x500	€ 140,00
☐	93056	Varilite Stratus	500x450	€ 140,00
☐	93057	Varilite Stratus	500x500	€ 140,00


 Varilite Zoid™
PSV

Low-profile, adjustable skin protection performance
cushion for active wheelchair users

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)	PRICE
☐	93031	Varilite Zoid™ PSV	300x300	€ 250,00
☐	93032	Varilite Zoid™ PSV	300x350	€ 250,00
☐	93033	Varilite Zoid™ PSV	350x350	€ 250,00
☐	93034	Varilite Zoid™ PSV	350x400	€ 250,00
☐	93100	Varilite Zoid™ PSV	350x450	€ 250,00
☐	93101	Varilite Zoid™ PSV	380x380	€ 250,00
☐	93102	Varilite Zoid™ PSV	380x430	€ 250,00
☐	93035	Varilite Zoid™ PSV	400x400	€ 250,00
☐	93036	Varilite Zoid™ PSV	400x450	€ 250,00
☐	93103	Varilite Zoid™ PSV	400x500	€ 250,00
☐	93104	Varilite Zoid™ PSV	430x430	€ 250,00
☐	93105	Varilite Zoid™ PSV	450x400	€ 250,00
☐	93037	Varilite Zoid™ PSV	450x450	€ 250,00
☐	93106	Varilite Zoid™ PSV	450x500	€ 250,00
☐	93107	Varilite Zoid™ PSV	500x400	€ 250,00
☐	93108	Varilite Zoid™ PSV	500x450	€ 250,00
☐	93109	Varilite Zoid™ PSV	500x500	€ 250,00



Varilite Junior

Ideal for children with a high degree
of tissue breakdown and symmetric
positioning needs

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)	PRICE
☐	93588	Varilite Junior	250x250	€ 250,00
☐	93599	Varilite Junior	275x325	€ 250,00
☐	93600	Varilite Junior	300x300	€ 250,00

5 Backrest

Backrest bracket

	CODE	DESCRIPTION	WIDTH x DEPTH (mm)	PRICE
☐	80298	Backrest bracket adjustable not in combination armrest 97471	seat width adjustable in steps 380 - 400 - 450 - 500 - 550 mm	€ 180,00
☐	80355	Backrest bracket adjustable Extra long tubes for mounting other back systems not in combination armrest 97471	seat width adjustable in steps 380, 400, 450, 500, 550 mm	€ 230,00
☐	80298-1	Backrest bracket adjustable JUNIOR not in combination armrest 97471	seat width adjustable in steps 250 - 300 - 350 mm	€ 260,00
☐	97485	Backrest Combi low (NEW), Compatible with Varilite Icon Back	seat width adjustable in steps 380 - 400 - 450 - 500 - 550 mm	410mm € 180,00
☐	97486	Backrest Combi Medium (NEW), Compatible with Varilite Icon or Netti back	seat width adjustable in steps 380 - 400 - 450 - 500 mm	51mm € 230,00
☐	97499	Backrest Combi Medium Large (NEW), Compatible w.	seat width adjustable in steps	51mm € 230,00
☐	97487	Backrest Combi High (NEW), Compatible with Varilite Icon or Netti back	seat width adjustable in steps 380 - 400 - 450 - 500 mm	61mm € 240,00
☐	97500	Backrest Combi High Large (NEW), Compatible with Varilite Icon or Netti back	seat width adjustable in steps 450 - 500 - 550 mm	61mm € 240,00



80298


 80298-1
Junior

 97485
Only in combination
with
Varilite Icon Back

 97486 - 97499
Compatible with
Varilite Icon or Netti
back

 97487 - 97500
Compatible with
Varilite Icon or Netti
back

Backrest system Netti (only in combination w. combi Medium or Combi High back)



97474
 Netti back straps
 To combine w.
 Combi Medium or
 Combi High back



97488
 Netti back raiser
 Back raiser to be used
 in combination with
 Combi High back



97498-1
 Netti extra strap
 For additional
 support in the lower
 back

CODE	DESCRIPTION	PRICE
☐ 97474	Standard Strap for backrest Mid. Used with 97486 and 97499	€ 80,00
☐ 97474	Standard Strap backrest Mid	
☐ 97488	incl. strap back top Used with 97487 and 97500	€ 100,00
☐ 97498-1	Strap backrest low, incl. brackets Used if an extra single strap is needed	€ 30,00



Netti Uno Back
 Good and firm back
 cushion with lumbar
 support. Steplessly
 adjustable in height



Netti Smart
 For users with
 higher requirements
 for comfortMultiple
 foam layers
 Steplessly adjustable
 in height



Netti Stabil
 For users with higher
 requirements for
 comfort
 Good lumbar support.
 Works well in
 combination with the
 Stabil lateral support.
 Multiple foam layers
 allows faster sink-in
 Steplessly adjustable in
 height



Netti Super Stabil
 For more comfort
 and stability
 Good lumbar
 support. Works well
 in combination with
 the Stabil lateral
 support
 Comes with two
 pads for adjustment
 of lumbar support
 and side support as
 well as adaptation
 for light scoliosis.



Netti Kyphotic
 For users with
 curved spine
 Support along the
 spine
 Multiple foam layers
 - increased sink-in
 Low lateral supports.
 Can be individually
 fitted steplessly in
 height

☐	N344702	Netti Uno Back 3D Mid	400	Multiple foam layers	€ 195,00
☐	N344703	Netti Uno Back 3D Mid	450		€ 195,00
☐	N344704	Netti Uno Back 3D Mid	500		€ 195,00
☐	N344705	Netti Uno Back 3D Mid	550		€ 195,00
☐	N344902	Netti Uno Back Easy Care Mid	400		€ 195,00
☐	N344903	Netti Uno Back Easy Care Mid	450		€ 195,00
☐	N344904	Netti Uno Back Easy Care Mid	500		€ 195,00
☐	N344905	Netti Uno Back Easy Care Mid	550		€ 195,00
☐	N357305	Netti Smart 3D Mid	380		€ 195,00
☐	N357306	Netti Smart 3D Mid	400		€ 195,00
☐	N357307	Netti Smart 3D Mid	430		€ 195,00
☐	N357308	Netti Smart 3D Mid	450		€ 195,00
☐	N357309	Netti Smart 3D Mid	500		€ 195,00
☐	N357310	Netti Smart 3D Mid	550		€ 305,00
☐	N358305	Netti Smart Easy Care Mid	380		€ 195,00
☐	N358306	Netti Smart Easy Care Mid	400		€ 195,00
☐	N358307	Netti Smart Easy Care Mid	430		€ 195,00
☐	N358308	Netti Smart Easy Care Mid	450		€ 195,00
☐	N358309	Netti Smart Easy Care Mid	500		€ 195,00
☐	N358310	Netti Smart Easy Care Mid	550		€ 305,00
☐	N357321	Netti Smart 3D Tall	380		€ 228,00
☐	N357322	Netti Smart 3D Tall	400		€ 228,00
☐	N357323	Netti Smart 3D Tall	430		€ 228,00
☐	N357324	Netti Smart 3D Tall	450		€ 228,00
☐	N357325	Netti Smart 3D Tall	500		€ 228,00
☐	N357328	Netti Smart one-piece foam 3D Tall	550		€ 367,00
☐	N358321	Netti Smart Easy Care Tall	380		€ 228,00
☐	N358322	Netti Smart Easy Care Tall	400		€ 228,00
☐	N358323	Netti Smart Easy Care Tall	430		€ 228,00
☐	N358324	Netti Smart Easy Care Tall	450		€ 228,00
☐	N358325	Netti Smart Easy Care Tall	500		€ 228,00
☐	N358328	Netti Smart one-piece foam Easy Care Tall	550		€ 367,00
☐	N367305	Netti Stabil 3D Mid	380		€ 326,00

☐	N367306	Netti Stabil 3D Mid	400	€ 326,00
☐	N367307	Netti Stabil 3D Mid	430	€ 326,00
☐	N367308	Netti Stabil 3D Mid	450	€ 326,00
☐	N367309	Netti Stabil 3D Mid	500	€ 326,00
☐	N367310	Netti Stabil 3D Mid	550	€ 424,00
☐	N368305	Netti Stabil Easy Care Mid	380	€ 326,00
☐	N368306	Netti Stabil Easy Care Mid	400	€ 342,00
☐	N368307	Netti Stabil Easy Care Mid	430	€ 342,00
☐	N368308	Netti Stabil Easy Care Mid	450	€ 342,00
☐	N368309	Netti Stabil Easy Care Mid	500	€ 342,00
☐	N368310	Netti Stabil Easy Care Mid	550	€ 424,00
☐	N367321	Netti Stabil 3D Tall	380	€ 362,00
☐	N367322	Netti Stabil 3D Tall	400	€ 362,00
☐	N367323	Netti Stabil 3D Tall	430	€ 362,00
☐	N367324	Netti Stabil 3D Tall	450	€ 362,00
☐	N367325	Netti Stabil 3D Tall	500	€ 362,00
☐	N367331	Netti Stabil 3D Tall	550	€ 411,00
☐	N368321	Netti Stabil Easy Care Tall	380	€ 362,00
☐	N368322	Netti Stabil Easy Care Tall	400	€ 362,00
☐	N368323	Netti Stabil Easy Care Tall	430	€ 362,00
☐	N368324	Netti Stabil Easy Care Tall	450	€ 362,00
☐	N368325	Netti Stabil Easy Care Tall	500	€ 362,00
☐	N368328	Netti Stabil one-piece Easy Care Tall	550	€ 436,00
☐	N377302	Netti Super Stabil 3D Mid	380	€ 368,00
☐	N377303	Netti Super Stabil 3D Mid	400	€ 368,00
☐	N377304	Netti Super Stabil 3D Mid	430	€ 368,00
☐	N377305	Netti Super Stabil 3D Mid	450	€ 368,00
☐	N377306	Netti Super Stabil 3D Mid	500	€ 368,00
☐	N377307	Netti Super Stabil 3D Mid	550	€ 505,00
☐	N378302	Netti Super Stabil Easy Care Mid	380	€ 387,00
☐	N378303	Netti Super Stabil Easy Care Mid	400	€ 387,00
☐	N378304	Netti Super Stabil Easy Care Mid	430	€ 387,00
☐	N378305	Netti Super Stabil Easy Care Mid	450	€ 387,00
☐	N378306	Netti Super Stabil Easy Care Mid	500	€ 387,00
☐	N378307	Netti Super Stabil Easy Care Mid	550	€ 505,00
☐	N377321	Netti Super Stabil 3D Tall	380	€ 424,00
☐	N377322	Netti Super Stabil 3D Tall	400	€ 424,00
☐	N377323	Netti Super Stabil 3D Tall	430	€ 424,00
☐	N377324	Netti Super Stabil 3D Tall	450	€ 424,00
☐	N377325	Netti Super Stabil 3D Tall	500	€ 424,00
☐	N377331	Netti SuperStabil one-piece foam 3D Tall	550	€ 535,00
☐	N378321	Netti Super Stabil Easy Care Tall	380	€ 424,00
☐	N378322	Netti Super Stabil Easy Care Tall	400	€ 424,00
☐	N378323	Netti Super Stabil Easy Care Tall	430	€ 424,00
☐	N378324	Netti Super Stabil Easy Care Tall	450	€ 424,00
☐	N378325	Netti Super Stabil Easy Care Tall	500	€ 424,00
☐	N378331	Netti SuperStabil one-piece foam Easy Care Tall	550	€ 535,00
☐	N206301	Netti kyphotic 3D/Mars Mid	350	€ 436,00
☐	N206302	Netti kyphotic 3D/Mars Mid	400	€ 436,00
☐	N206303	Netti kyphotic 3D/Mars Mid	450	€ 436,00
☐	N206304	Netti kyphotic 3D/Mars Mid	500	€ 436,00

Backrest system Varilite Icon™ (only in combination w. Combi Low back)



Icon™ Back System
Low

Combines innovative
Varilite™ hardware
with a sleek, low
profile look



Icon™ Back System
Mid

Provides support
and positioning to
improve trunk
control and
decrease fatigue



Icon™ Back System
Tall

Well suited for users
needing additional
support or with a tilt-
recline system



Icon™ Back System
Deep

Ideal for users
needing maximum
support

☐	93560	Varilite Icon™ Low, incl. adapter set 93567	350	€ 360,00
☐	93561	Varilite Icon™ Low, incl. adapter set 93567	380	€ 360,00
☐	93562	Varilite Icon™ Low, incl. adapter set 93567	400	€ 360,00
☐	93563	Varilite Icon™ Low, incl. adapter set 93567	430	€ 360,00
☐	93564	Varilite Icon™ Low, incl. adapter set 93567	450	€ 360,00
☐	93565	Varilite Icon™ Low, incl. adapter set 93567	500	€ 360,00
☐	93587	Varilite Icon™ Mid, incl. adapter set 93567	250	€ 420,00
☐	93544	Varilite Icon™ Mid, incl. adapter set 93567	300	€ 360,00

☐	93545	Varilite Icon™ Mid, incl. adapter set 93567	350	€ 360,00
☐	93546	Varilite Icon™ Mid, incl. adapter set 93567	380	€ 360,00
☐	93547	Varilite Icon™ Mid, incl. adapter set 93567	400	€ 360,00
☐	93548	Varilite Icon™ Mid, incl. adapter set 93567	430	€ 360,00
☐	93549	Varilite Icon™ Mid, incl. adapter set 93567	450	€ 360,00
☐	93550	Varilite Icon™ Mid, incl. adapter set 93567	500	€ 360,00
☐	93551	Varilite Icon™ Mid, incl. adapter set 93567	550	€ 420,00
☐	93552	Varilite Icon™ Mid, incl. adapter set 93567	600	€ 420,00
☐	93553	Varilite Icon™ Tall, incl. adapter set 93567	300	€ 390,00
☐	93554	Varilite Icon™ Tall, incl. adapter set 93567	350	€ 390,00
☐	93555	Varilite Icon™ Tall, incl. adapter set 93567	380	€ 390,00
☐	93556	Varilite Icon™ Tall, incl. adapter set 93567	400	€ 390,00
☐	93557	Varilite Icon™ Tall, incl. adapter set 93567	430	€ 390,00
☐	93558	Varilite Icon™ Tall, incl. adapter set 93567	450	€ 390,00
☐	93559	Varilite Icon™ Tall, incl. adapter set 93567	500	€ 390,00
☐	93586	Varilite Icon™ Deep, incl. adapter set 93567	250	€ 500,00
☐	93538	Varilite Icon™ Deep, incl. adapter set 93567	350	€ 500,00
☐	93539	Varilite Icon™ Deep, incl. adapter set 93567	380	€ 500,00
☐	93540	Varilite Icon™ Deep, incl. adapter set 93567	400	€ 500,00
☐	93541	Varilite Icon™ Deep, incl. adapter set 93567	430	€ 500,00
☐	93542	Varilite Icon™ Deep, incl. adapter set 93567	450	€ 500,00
☐	93543	Varilite Icon™ Deep, incl. adapter set 93567	500	€ 500,00

Backrest adjustments and options

CODE	DESCRIPTION	adjustment range	PRICE
☐	80363	Mechanical backrest adjustment, quick release	€ 110,00
☐	80577	Mechanical backrest adjustment, gas spring	€ 52,00
☐	97235	Electric backrest adjustment (100mm actuator)	€ 120,00
☐	81137	Electric backrest adjustment (150mm actuator)	€ 280,00



80363



80577



97235

Lateral support

CODE		DESCRIPTION					PRICE
for seat width							
30cm 93254 	35-38cm 93255 	40-43cm 93256 	45cm 93257 	50-60cm 93258 	Varilite lateral support fixed, set only in combination with Varilite Icon™	Choose underneath Cushion for thigh/lateral support for this. € 97,00	
30cm 93239 	35-38cm 93240 	40-43cm 93241 	45cm 93242 	50-60cm 93243 	Varilite lateral support swing away, set only in combination with Varilite Icon™	short Choose underneath Cushion for thigh/lateral support for this. € 250,00	
30cm 93249 	35-38cm 93250 	40-43cm 93251 	45cm 93252 	50-60cm 93253 	Varilite lateral support swing away, set only in combination with Varilite Icon™	High Choose underneath Cushion for thigh/lateral support for this. € 250,00	
30cm 93244 	35-38cm 93245 	40-43cm 93246 	45cm 93247 	50-60cm 93248 	Varilite lateral support swing away, set only in combination with Varilite Icon™	long Choose underneath Cushion for thigh/lateral support for this. € 250,00	
	93262	Cushion for thigh/lateral support, pcs 7,5x10 cm, set				€ 50,00	
	93263	Cushion for thigh/lateral support, pcs 10x13 cm, set				€ 50,00	
	93264	Cushion for thigh/lateral support, pcs 13x15 cm, set				€ 50,00	



Fixed



Swing away short



Swing away High



Swing away long


Head support

CODE	DESCRIPTION	Size (mm) / Shape	Cover	PRICE
	80337 TA Head support incl. bracket (80288)	225 x 125		€ 114,00
	80829 TA Head support High incl. bracket (80288-2)	175 x 95		€ 130,00
	80267 TA Head support with support pad			€ 210,00



80267

	M1086677 Meyra Multifunctional 3D head support, including bracket	Only for combi back rest mid + high vers.	textile black	€ 350,00
	M1086678 Meyra Multifunctional 3D head support, including bracket		Leather black	€ 350,00



M1086677
M1086678

	N42370	Netti head support A, incl. bracket (97518 & 80047)	curved	3D	€ 165,00
	N42369	Netti head support A, incl. bracket (97518 & 80047)	curved	Easy Care	€ 115,00
	N87398	Netti head support A with forehead strap, incl. bracket (97518 & 80047)	curved	3D	€ 178,00
	N87397	Netti head support A with forehead strap, incl. bracket (97518 & 80047)	curved	Easy Care	€ 178,00
	N44579	Netti head support A Mini, incl. bracket (97518 & 80047)	150 x 200	3D	€ 203,00
	N42361	Netti head support B, incl. bracket (97518 & 80047)	small	3D	€ 115,00
	N42360	Netti head support B, incl. bracket (97518 & 80047)	small	Easy Care	€ 115,00
	N42352	Netti head support C, incl. bracket (97518 & 80047)	large	3D	€ 115,00
	N42351	Netti head support C, incl. bracket (97518 & 80047)	large	Easy Care	€ 115,00
	N81021	Netti head support D, incl. bracket (97518 & 80047)	comfort	3D	€ 143,00
	N82983	Netti head support D, incl. bracket (97518 & 80047)	comfort	Easy Care	€ 143,00
	N94233	Netti head support E, incl. bracket (97518 & 80047)	adjustable sides	3D	€ 216,00
	N94032	Netti head support E, incl. bracket (97518 & 80047)	adjustable sides	Easy Care	€ 216,00
	N87414	Netti head support F, incl. bracket (97518 & 80047)	with chin support	3D	€ 216,00
	N87413	Netti head support F, incl. bracket (97518 & 80047)	with chin support	Easy Care	€ 216,00



Netti head support
A



Netti support A with
forehead strap



Netti head support
Mini



Netti head support B



Netti head support C



Netti head support
D




Netti head support
E



Netti head support F

☐	N45365	Netti Savant Pediatric head support	120 mm	€ 469,00
☐	N45366	Netti Savant Small head support	130 mm	€ 469,00
☐	N45367	Netti Savant Adult head support	150 mm	€ 469,00
☐	N45369	Horizontal bar for Netti Savant head support	150 mm	€ 126,00
☐	N45370	Horizontal bar for Netti Savant head support	250 mm	€ 126,00
☐	N45368	Netti Axion head support adapter 60° to combine w. Savant head support		€ 514,00
☐	N45378	Netti Axion head support adapter 90° to combine w. Savant head support		€ 543,00
☐	N45373	Netti Savant side Pad		€ 48,00
☐	N45398	Netti Savant forehead band		€ 126,00



Netti Savant head support



Netti Axion head support
adapter



Netti Savant side pad

Thigh support

	CODE	DESCRIPTION	PRICE
☐	93262 80328	Thigh support, incl. bracket, cushion 75 x 100 mm, set	set price € 270,00
☐	93263 80328	Thigh support, incl. bracket, cushion 100 x 130 mm, set	set price € 270,00
☐	93264 80328	Thigh support, incl. bracket, cushion 130 x 150 mm, set	set price € 270,00

Belts & support

CODE	DESCRIPTION				PRICE
☐ 93231	Chest belt				€ 27,00
	XS	S	M	L	
☐ 93221	☐	☐	☐	☐	2-point hip belt € 46,00
☐ 93225	☐	☐	☐	☐	4-point hip belt € 58,00
☐ 93235	☐	☐	☐	☐	Band-It, Ankle, set € 74,00
☐ 93272	☐	☐	☐	☐	Contoured chest harness € 83,00
☐ 80765	Bracket for chest harness				€ 45,00
☐ 97157	Belt, automatic roll-up				€ 200,00

6 Arm support

CODE	DESCRIPTION		PRICE
☐ 97471	Arm support flip-up with click function, set	height adjustable from 230 to 320 mm	€ 325,00
☐ 80289-2	Arm support flip-up with front support lockable, set	height adjustable from 185 to 285 mm	€ 250,00
☐ 80302	Arm support flip-up, set	height adjustable from 185 to 285 mm	€ 250,00
☐ 80743	Stop for Flip up arm support pcs.in comb. w. 80302 & 80289-2	prevents arm support from going further back	€ 7,00
☐ 80347	Arm support bracket rigid, set	adjustable from 100 to 280 mm from the seat plate	€ 170,00



80289-2



80302



80347



97471

Options for armsupport

CODE	DESCRIPTION		PRICE
☐ 80320	Arm pad cushion, right, short 320mm	length 320mm	€ 25,00
☐ 80335	Arm pad cushion, right, long 430mm	length 430mm	€ 28,00
☐ 80336	Arm pad cushion, left, short 320mm	length 320mm	€ 25,00
☐ 80321	Arm pad cushion, left, long 430mm	length 430mm	€ 28,00
☐ 80651	Aluminium edge for arm support cushion short		€ 23,00
☐ 80652	Aluminium edge for arm support cushion long		€ 26,00
☐ 80830	Swing bracket for arm pad right	for armrest flip-up 80302	€ 39,00
☐ 80830-1	Swing bracket for armpad left	for armrest flip-up 80302	€ 39,00
☐ 80994	Purse -leather bag for arm support		€ 140,00
☐ 20248	Purse -pocket for arm support		€ 67,00



80320/80336



80335/80321



80651/80652



80830/80830-1



80994



20248

7 Leg support and Footplate

Two-piece leg support

CODE	DESCRIPTION	PRICE
☐ 80700-1	Manual elevating leg support, right	€ 280,00
☐ 80701-1	Manual elevating leg support, left	€ 280,00
☐ 80628	Electric elevating leg support, right	€ 450,00
☐ 80629	Electric elevating leg support, left	€ 450,00
☐ 80716	Amputation elevating leg support right	€ 170,00
☐ 80717	Amputation elevating leg support left	€ 170,00
☐ 80721	Knee pad standard right	for electric separate leg support, right 80628 € 96,00
☐ 80730	Knee pad standard left	for electric separate leg support, left 80629 € 96,00
☐ 80631-3	Lower tube for leg support straight + 100mm	€ 20,00
☐ 80631-2	Lower tube for leg support curved + 100mm	€ 24,00
☐ M293	Fixed angle leg support, set	incl. bracket adaption for MeyLife legrest (97449) € 395,00
☐ M4092	Manual elevating leg support, set incl. bracket adaption for MeyLife legrest	incl. bracket adaption for MeyLife legrest (97449) € 1.035,00
☐ M2086	Electric elevating leg support, set incl. bracket adaption for MeyLife legrest	incl. bracket adaption for MeyLife legrest (97449) € 2.033,90



80700-1/80701-1



80628/80629



80716/80717



80721 / 80730



M293

Central leg support with single footplate

CODE	DESCRIPTION	PRICE
☐ 95083	Center manual leg support adjustable	€ 300,00
☐ 80357	Center power leg support with calf pads	mechanically length adjustment (via screws) 340-430mm + 150mm electrical length adjustment (via controller) € 820,00
☐ 81124	Center power leg support with calf pads for separate footplates	mechanically length adjustment (via screws) 340-430mm + 150mm electrical length adjustment (via controller) € 820,00
☐ 81128 81017	Center power leg support with calf pads incl. electric flip-up footplate	mechanically length adjustment (via screws) 340-430mm + 150mm electrical length adjustment (via controller) € 3.000,00
☐ 80571	Center power leg support short with calfpads	mechanically length adjustment (via screws) 250-310mm + 80mm electrical length adjustment (via controller) € 820,00
☐ 81126	Center power leg support short with calfpads for separate footplates	mechanically length adjustment (via screws) 250-310mm + 80mm electrical length adjustment (via controller) € 820,00



95083



81017
Footplate flip-up
electric



80357

Footplate for central leg support

CODE	DESCRIPTION	PRICE
☐ 80373-2	Footplate moulded standard complete with edge. Not good with a FWD	€ 0,00
☐ 80539-1	Footplate with cutout moulded complete with edge	€ 0,00
☐ 80818	Footplates divided	€ 520,00



80373-2



80539-1



80818

Options for central leg support

CODE	DESCRIPTION	PRICE
☐ 80690	Knee protection for center leg support	€ 46,00



80690

8 Lighting

CODE	DESCRIPTION	PRICE
☐ 97223	Active LED lighting	€ 140,00
☐ ---	Without lighting	€ 0,00

9 Control module and accessories

CODE	DESCRIPTION	PRICE
☐ 95174	Controller, JSM, LED Entry level	€ 210,00
☐ 95253	Controller, JSM-L LED	€ 240,00
☐ 95172-1	Controller, CJSM, SW (direct actuator)	€ 320,00
☐ 95252-1	Controller, CJSM-L-SW(direct actuator)	€ 320,00
☐ 95501	Controller CJSM2-L	€ 480,00
☐ 95524	Controller, CJSM2, Bluetooth	€ 590,00
In combination with CJSM or CJSM 2, choose your preferred language		
☐	Danish	
☐	English	
☐	German	
☐	French	



95174/95253
JSM LED



95172-1/95252-1
CJSM



95501/95524
CJSM-2

Control module bracket

CODE	DESCRIPTION	PRICE
☐ 80407	Right swingaway bracket for control module right	€ 76,00
☐ 80408	Left swingaway bracket for control module left	€ 76,00
☐ 80838	Drop-Down bracket universal	€ 120,00



80407/80408



80838

Control module knob

	CODE	DESCRIPTION	PRICE
☐	95153	Standard	€ 0,00
☐	95154	Mushroom M48	€ 41,00
☐	95155	Ball 38 mm	€ 28,00
☐	95299	Ball 32mm	€ 26,00
☐	95306	Chin C38	€ 26,00
☐	95298	Cone S30	€ 11,00
☐	95305	Stick S80	€ 41,00
☐	95301	T-bar	€ 50,00
☐	95300	Softball 70mm	€ 41,00
☐	Select color	Golfball	€ 10,00



95153



95154



95155/95299



95306



95298



95305



95301



95300



White
80653



Orange
80653-4



Pink
80653-2



Yellow
80653-3



Football
80653-7



Football
80653-8

Special operating modules

	CODE	DESCRIPTION	PRICE
☐	97342	TA hand control for CJS, incl. bracket (80678)	€ 98,00
☐	97343	TA hand control for CJS2, incl. bracket (80678)	€ 98,00
☐	95213-1	Attendant controller	€ 520,00
☐	M1087331	Attendant controller for Combi back complete, left	€ 600,00
☐	M1087330	Attendant controller for Combi back complete, right	€ 600,00
☐	80903	1 hand Attendant bracket swing	without controller 95213-1 € 200,00
☐	80517	1 hand attendant bracket IW / iQ	without controller 95213-1 € 68,00
☐	80533	TA 2-hand attendant control, with emergency stop	without controller 95213-1 € 660,00



97342/97343

10 Accessories

	CODE	DESCRIPTION	PRICE
☐	80163	Tray, adult, incl. bracket (80518)	€ 303,00
☐	80164	Tray, junior, incl. bracket (80518)	€ 303,00
☐		Half tray, on request and customized	
☐	80175	Bag hooks	€ 33,00



80162/80164



80175