

HURRICANE PRO

1.880 CUSTOM DESIGN



Permissible user weight
100 kg



MEYRA®



Adjustable oversize axle tube for ideal gravity setting



Footplate individually adjustable



Double support castor wheels standard



Aluminium side panels screw-retained

The precision device

EFFICIENCY

- An economic acquisition for institutions offering wheelchair sports for children and youths, and for clubs seeking to support young players
- Multiple use
- Fully equipped, incl. double support castors
- Optimum cost-performance and quality ratio

THERAPEUTIC BENEFIT

- Individual adjustments can be carried out without great effort, as the relevant key data are already defined in the framework concept
- The seating position and an optimal centre of gravity provide stability during mobility through ergonomics and little physical effort for locomotion
- Outstanding handling, optimum smoothness of ride, precise straight-line stability and excellent manoeuvrability promote sports activity in both rehabilitation and leisure
- Adjustability of the centre of gravity and seat width

USE

- For sports clubs, schools and other institutions dedicated to wheelchair sports
- For users who play wheelchair basketball as a grassroots or leisure sport
- For wheelchair basketball newcomers

ORDER / CUSTOM DESIGN

COST ESTIMATE

Please fax the completed order form to: +49 5733 922 9311 Customer service: +49 5733 922 311

Billing address / Customer ref. no.:
Delivery address:
Consignment / Comments:

MEYRA®



 Permissible user weight
100 kg

HURRICANE PRO 1.880 CUSTOM DESIGN

All prices subject to VAT at the prevailing rate! Valid from 01.01.2023 – 31.12.2023.
Our Terms and Conditions as stated at www.meyra.de/AGB shall apply.
Separate delivery terms for the HURRICANE 1.880 apply in addition to the normal terms and conditions.

00 Basic equipment

MODEL	NAME	
<input type="checkbox"/> 1880	HURRICANE PRO Basketball wheelchair	Drive wheels with Jumbo rim, Speed Air tyres, stainless steel hand rim, double support castors, aluminium side panels screw-retained, adjustable footplate, round deflector bracket, adjustable oversize axle tube, adjustable back strap, seat cushion, strapping belt.

Please select by ticking the relevant boxes

01 Frame version

CODE	NAME
<input type="checkbox"/> 353	24"
<input type="checkbox"/> 354	26"

Equipment for frame version

NAME	CODE 353 (24") INCLUDES	CODE 354 (26") INCLUDES
Frame length	650 mm	700 mm
Specify back height when ordering, otherwise maximum measurement	290 – 330 mm	260 – 300 mm
Front seat height	480 mm	530 mm
Rear seat height	380 mm	480 mm
Axle position, infinitely variable adjustment	Adjustment range from 60 to 170 mm (specify when ordering, otherwise maximum measurement)	Adjustment range from 60 to 170 mm (specify when ordering, otherwise maximum measurement)
Lower leg length, measured without cushion	310 mm – 430 mm	360 mm – 480 mm
Backrest angle, fixed	90° to road surface	90° to road surface
Side panel aluminium	Height: 237 mm Length: 405 mm	Height: 164 mm Length: 405 mm
Total width	Seat width + 440 mm 24"	Seat width + 460 mm 26"
Drive wheel size	Quick release axle 12.7 mm Aluminium hub Jumbo rim	Quick release axle 12.7 mm Aluminium hub Jumbo rim
Hand rim	V2A, 6 brackets	V2A, 6 brackets
Tyres, Schwalbe	Speed Air	Speed Air
Wheel camber, fixed	16°	16°
Footboard	height and angle-adjustable	height and angle-adjustable
Skater wheel	80 mm	80 mm
Back strap	adaptable	adaptable
Seat belt	yes	yes
Strapping belt HURRICANE The Original	yes	yes
Double support castors Version, welded	60 mm skater	60 mm skater
Deflector round, welded	yes	yes

02 Seat width

CODE	NAME
<input type="checkbox"/> 32	320 mm
<input type="checkbox"/> 36	360 mm
<input type="checkbox"/> 40	400 mm

Please select by ticking the relevant boxes

03 Seat depth

CODE	NAME
<input type="checkbox"/> 392	320 mm
<input type="checkbox"/> 366	360 mm
<input type="checkbox"/> 615	400 mm

04 Seat cushion

CODE	NAME
<input type="checkbox"/> 740	Seat cushion 50 mm soft
<input type="checkbox"/> 720	Seat cushion 50 mm hard
<input type="checkbox"/> 726	Additional cushion 50 mm extra hard (Natio hardness)
<input type="checkbox"/> 740	Additional cushion 50 mm soft
<input type="checkbox"/> 720	Additional cushion 50 mm hard
<input type="checkbox"/> 996	Wedge cushion 50 mm to 20 mm (for changing tilt)

05 Strapping belt

CODE	NAME	SIZE XS (circumference 800 – 900 mm)	SIZE S (circumference 900 – 1,000 mm)	SIZE M (circumference 1,000 – 1,100 mm)	SIZE L (circumference 1,100 – 1,200 mm)	SIZE XL (circumference 1,200 – 1,300 mm)
<input type="checkbox"/> 362	Strapping belt	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> 362	Additional Strapping belt	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
<input type="checkbox"/> 361	Strapping belt Body Fix (ratchet system)					

06 Optional equipment (mostly for self-assembly)

CODE	NAME
<input type="checkbox"/> 490	Hand protectors and spoke guards
<input type="checkbox"/> 597	Leg protection pad
<input type="checkbox"/> 316	Leg fixation (for mounting on front frame)
<input type="checkbox"/> 445	Titanium hand rim, 6 brackets
<input type="checkbox"/> 4953	Spinergy wheel (one pair, surcharge for new vehicle)

Date:	Signature:	As at: 2023-01-01
-------	------------	-------------------