

EUROCHAIR² HD / INDIVI- DUAL




Authorised user weight
now with limit of 300 kg



MEYRA[®]



Highly durable



One-piece footrest



Double cross offers more stability



Reinforced seat belt

Maximum stability

COST EFFECTIVENESS

- as a fully equipped chair for active, long-term use in everyday life, this model is versatile and adjustable to the individual user
- the specific use of rigid connecting components in the frame gives the folding wheelchair the stability it needs for this purpose
- long-term supply of spare parts thanks to extensive network of specialist retailers
- ideal for reuse thanks to simple adjustment

THERAPEUTIC BENEFITS

- easier to use thanks to numerous adjustment settings for the chassis the grip angle is ergonomically designed to relieve the shoulders
- improved effectiveness increases independence and facilitates optimum utilisation of existing propulsion power, easier movement, high manoeuvrability
- a horizontal push bar and/or one-piece footrest ensure that the user receives maximum support while sitting

APPLICATIONS

- particularly suitable for heavy persons
- designed for long-term care with high adjustability for independent mobility
- predominantly for domestic use, but also useful in rehabilitation clinics for achieving therapeutic goals
- can be used with various forms of propulsion

SPECIAL BUILD

QUOTATION

Please fax your completed order form to: +49 (0) 5733 922 9311 Customer service: +49 (0) 5733 922 311

Invoice address/customer no.:

Delivery address:

Consignment/remarks:

Address: Ms Mr

Last Name:

First Name:

Street/Number.:

Postal CODE/City:

Date of birth:

Height: _____ Weight: _____



Authorised user weight up to 300 kg

EUROCHAIR² HD / INDIVIDUAL

Reinforced design for user weight of up to 300 kg

- Reinforced rear strap
 - Reinforced seat belt for SW 580 and 650
 - Reinforced rear side guard intake
- Only in combination with the following features:
- CODE 54 one-piece footrest
- Wheelchair not suitable for transporting persons.



MEYRA

Example model

All prices exclude VAT at the statutory rate. Valid from 01.01.2019 to 31.12.2019.
Our general terms and conditions apply, which can be found on www.meyra.de/en/gtc.



Technical data

Seat width	700	750	800
Width when ready for use	920	970	1,020
Width folded	310		
Seat depth	460 / 500		
Seat height front/rear	500 / 475 // 475 / 450		
Overall height	960		
Back height	430		
Armrest height	230 – 330		
Length with footrests	max. 1,140		
Length without footrests	max. 880		

User weight	300
Empty weight	34
Permitted overall weight	334
Transport weight	27,5
Rear wheel	610 (24")

Dimensions in mm (1 inch = 25,4 mm), Weights in kg (1 pd = 0,45 kg), Speed in km/h (1 mph = 1,61 km/h), Tolerance ± 10 mm / 0,39 inch
The specifications are subject to constructional changes.

Caution!
Please note when using public transport (e.g. trains): the dimensions of the wheelchair may exceed the transport company guidelines. Please enquire beforehand.

Please tick your choice in the boxes

Standard version: In case you do not make another selection we will generally deliver the standard version.

01 Seat width (SW)

	CODE	DESCRIPTION	
<input type="checkbox"/>	3368	SW 700 mm	User weight up to 300 kg
<input type="checkbox"/>	3370	SW 750 mm	User weight up to 300 kg
<input type="checkbox"/>	3372	SW 800 mm	User weight up to 300 kg

02 Seat depth (SD)

	CODE	DESCRIPTION
<input type="checkbox"/>	619	SD 460 mm
<input type="checkbox"/>	792	SD 500 mm

03 Seat height front (FSH) combined with castor and rear wheel (without seat cushion)

Please select your desired seat height by choosing a possible combination of seat height and wheel size (only 24 x 1 3/8" rear wheel).

CODE	Ø CASTOR IN MM	FSH IN MM:	
		475	500
154	Ø 142 soft solid rubber	<input type="checkbox"/>	
157/900	Ø 200 flexel		<input type="checkbox"/>

When FSH is selected then the standard castor and matching castor fork are automatically configured.

The rear seat height follows automatically from the standard seat inclination of approx. 4°, which accounts for a difference of approx. 30 mm.

04 Seat cushion

	CODE	DESCRIPTION
<input type="checkbox"/>	561	Seat cushion 60 mm
<input type="checkbox"/>	714	incontinence cover for seat cushions, grey

05 Hand rims and hand rim options

	CODE	DESCRIPTION	Ø
<input type="checkbox"/>	493	Stainless steel hand rim	24"
<input type="checkbox"/>	692	Hand rim gap of 15 mm	24"
<input type="checkbox"/>	693	Hand rim gap of 25 mm (only in combination with CODE 493)	

Please tick your choice in the boxes

Standard version: In case you do not make another selection we will generally deliver the standard version.

06 Tyre

CODE	DESCRIPTION
<input type="checkbox"/> 904	Puncture-resistant Flexel tyres

07 Brakes

CODE	DESCRIPTION
<input type="checkbox"/> 680/706	Pressure brake for driver and drum brake for accompanying person
<input type="checkbox"/> 986	Brake lever extension, attachable, 1 pair

08 Back strap height

CODE	DESCRIPTION	HEIGHT
<input type="checkbox"/> 442	Back with reinforced design	430 mm

09 Backrest

CODE	DESCRIPTION
<input type="checkbox"/> 622	Backrest with lumbar curve

Backrest options

CODE	COMPATIBLE WITH CODE...
<input type="checkbox"/> 814	Horizontal push bar, stabilises back in the case of heavier loads

10 Side guard and armrest

Please note the marked functional features!

CODE	DESCRIPTION
<input type="checkbox"/> 106	Side guards height adjustable, armrest cushion depth adjustable

11 Leg rest



removable



swing-away



swing-up


CODE	DESCRIPTION
<input type="checkbox"/> 94	Without leg rests and footrests
<input type="checkbox"/> 93	Leg rests, removable and swing-away


Please tick your choice in the boxes


Standard version: In case you do not make another selection we will generally deliver the standard version.

12 Footrest

 flip-up

 height adjustable

 depth adjustable

 angle adjustable

CODE	DESCRIPTION	DEPTH	LSL
<input type="checkbox"/> 54	One-piece footrest, depth adjustable by 60 mm (in 4 steps), fold-up, height and angle adjustable	150	350 - 490

Required LSL in mm (measured without seat cushion)

Length adjustment by manufacturer of 350-520 mm available for calf depending on coding.

Footrest accessories

CODE	DESCRIPTION
<input type="checkbox"/> 822	Shoe positioning straps (pair)

13 Frame colour

CODE	DESCRIPTION
<input type="checkbox"/> 4919	Silverline

14 Seat and back strap colour

CODE	DESCRIPTION
<input type="checkbox"/> 287	Deep black, polyester

15 Wheelchair accessories

Hand and spoke guard

CODE	DESCRIPTION
<input type="checkbox"/> 489	Transparent

Anti-tip wheel

<input type="checkbox"/> 691	Anti-tip wheels, length adjustable, removable (1 pair)
------------------------------	--

Other accessories

<input type="checkbox"/> 585	Infusion stand "UNIVERSAL"
<input type="checkbox"/> 833	Safety belt with lock
<input type="checkbox"/> 819	Reflectors

Date:

Signature:

must notify us in writing of any complaints within one week of receipt of delivery. After expiry of this period, any liability for such defects is excluded. In the event of concealed defects, warranty claims shall only be valid if we are notified in writing within one week of discovery.

8.2. In cases of defective delivery, we shall have the right to choose to repair the defective item or to replace the defective item with a new delivery. The customer is obliged to provide us the defective product delivered free for inspection and repair. Should the repair or replacement delivery be unsuccessful, should they not be undertaken with a reasonable period or should we refuse to undertake them, the customer may withdraw from the contract or reduce the purchase price. A deadline shall not be required in those cases where it is not required by law.

8.3. The warranty shall be void if the customer modifies the delivery item without our consent or it is changed by a third party and it therefore becomes impossible or unreasonably difficult to rectify the defect. However, in each case, the customer must bear the additional costs resulting from the rectification of the defect.

8.4. The above regulations shall not affect any warranty that we assume on the part of the first user of the rehabilitation equipment.

8.5. Excluded from the warranty is the function-related wear of all items, modules, batteries and spare parts supplied by us as well as the inappropriate or improper storage, use or handling of the products.

8.6. Further claims by the customer, in particular for damages instead of performance and for reparation of any direct or indirect damage – including incidental or consequential damages, irrespective of the legal grounds – shall be excluded. This shall not be the case if MEYRA GmbH has concealed a legal or material defect maliciously or the damage is due to intent or gross negligence on the part of MEYRA GmbH, its legal representatives or vicarious agents or to any negligent breach of essential contractual obligations. Essential contractual obligations are understood to be those where their fulfilment enables the proper execution of the contract and on compliance with which the contracting parties regularly rely or may rely. In the case of damage to property and assets due to negligence, the liability of MEYRA GmbH shall however be restricted to an amount foreseeable in a typical contract. Furthermore, personal injury or damage to health caused by a breach of duty on the part of MEYRA GmbH, its legal representatives or vicarious agents are not excluded.

8.7. Liability under the *Produkthaftungsgesetz* [German Product Liability Act] is unaffected.

9. STATUTE OF LIMITATIONS

Claims by the customer relating to defects shall have a limitation period of two years from handover / delivery of the goods to the customer. With respect to replacement items and repairs, the limitation period shall be one year from handover / delivery of the goods to the customer. Excluded from this are claims for damages due to endangerment of life, personal injury or health and / or claims arising out of damages caused by gross negligence or intent on the part of MEYRA GmbH or its vicarious agents. In such instances, the statutory limitation periods shall apply.

9.1. REPLACEMENT PRODUCT SERVICE

Independent of the warranty, MEYRA GmbH offers customers the opportunity to exchange a defective product for a refurbished product. The following provisions shall apply to the products offered as replacements:

- Replacement products are refurbished and technically in order.
- The defective product must be returned free of charge to MEYRA GmbH within 15 working days.

The returned item shall become our property.

- Should the defective product not be returned, we shall invoice the customer for the replacement item at 95% of the sales price for a new product.
- The returned item must correspond to the replacement product in model and design. The product must also be reusable and may only show wear and tear commensurate with normal use.
- MEYRA GmbH shall determine whether the returned item fulfils the requirements specified above. Should these requirements not be fulfilled, we shall also issue the customer with an invoice for the amount specified above for the replacement product delivered less the residual value of the returned item.

10. PRODUCT SAFETY GUARANTEE

10.1. PRODUCT RECALL

If, due to defects in the goods delivered or to existing legal stipulations, it becomes necessary to have a product recalled, then the Client is obliged to cooperate in implementing such product recall in order to ensure that it can be conducted efficiently. In the event of a product recall, the Client is in particular obliged to inform MEYRA GmbH regarding the whereabouts of the delivery items in the delivery chain (from the Client to end consumer). To this end, the Client must continually document, using appropriate documentation methods, the progress of the goods and provide MEYRA GmbH with such documentation upon request. If the Client supplies the goods to resellers, then it must impose upon such resellers, via contractual provisions, the foregoing duties of cooperation intended to ensure the efficient implementation of a product recall.

10.2. The Client is obliged to have end customers instructed with regard to the handling and use of the delivered products and to have such instruction carried out by trained specialists. If the Client supplies the goods to resellers, then it must impose upon such resellers, via contractual provisions, the foregoing training obligations pertaining to the handling and use of the delivered goods.

11. LIABILITY FOR DESIGN CHANGES

It is important to note that more stringent legal provisions shall apply in the case of custom orders. Design changes to MEYRA GmbH products by the customer or a third party contracted by the customer shall only be permitted if they are in line with the safety requirements and our management has declared its consent in writing in advance. To this end, if requested, a modified model and design drawings are to be made available to us. Should design changes be undertaken without the prior written consent of our management and damage due to the changes is caused to a third party for which we are liable in the external relationship, the customer shall be required to release us in the internal relationship from claims by the third party.

12A. RETURN OF FINISHED PRODUCTS AND MODULES

Goods returned without an accompanying delivery note or invoice copy shall not be accepted. Credit equivalent to 80 % of the net value of the goods shall be granted for goods returned in their original packaging and in mint condition. Excluded from return are items delivered more than three months ago, custom orders, sanitary items, filled batteries and products with a net value of less than € 100.00. Also excluded are custom-made wheelchairs (e.g.: wheelchairs for children or adaptive wheelchairs). The transport risk shall be borne by the sender.

12B. ORDER CANCELLATIONS FOR FINISHED PRODUCTS

12B.1. Order cancellations without the prior agreement of MEYRA GmbH shall not be permitted.

12B.2. Cancellations shall be made in writing.

12B.3. Should an agreement be reached regarding the cancellation of finished products which are ready for delivery, the following shall apply:

- Rehabilitation products shall be credited in full.

In principle, orders for custom wheelchairs and sports wheelchairs cannot be cancelled.

- In the event of the cancellation of an order for adaptive wheelchairs, 20% of the net value shall be retained to cover costs.

The cancellation fee for electric wheelchairs and scooters is 5%. The cancellation fee for standard and lightweight wheelchairs is 10%.

12C. COLLECTION / DISPOSAL

Our prices do not include the costs for collection and disposal of complete old equipment from users other than private households. On request, and against reimbursement of the costs incurred, we are able to organise the collection and recycling / disposal of such equipment, provided it was purchased from us. Returns without the prior agreement of MEYRA GmbH are not possible.

13. USE OF PERSONAL DATA

We shall be entitled to save and process internally any personal data of the customer falling within the scope of the statutory provisions, in particular the General Data Protection regulation.

14. JURISDICTION, APPLICABLE LAW, SEVERABILITY CLAUSE

14.1. Unless expressly agreed to the contrary, the place of performance shall be the premises of MEYRA GmbH in Kalletal-Kalldorf.

14.2. The jurisdiction for all disputes arising out of the contractual relationship shall be determined by our headquarters in Kalletal-Kalldorf. We shall also be entitled to take action against the customer within its jurisdiction.

14.3. The law of the Federal Republic of Germany shall apply excluding the UN Convention on Contracts for the International Sale of Goods.

14.4. Should individual provisions of these General Terms and Conditions be wholly or partially void or ineffective, the validity of the remaining provisions shall remain unaffected.