iTRAVEL

1.054 Stock version

Allowed user weight

120 kg

MEYRA QUALITY TESTED



ORTOPEDIA

EINE MARKE VON MEYRA



Adjustable backrest



Comfort seat system



Flip-up armrests



Small pack size

Travel Wheelchair

ADAPTATION

- Weighing in at 22 kg, this is the lightweight among the electric wheelchairs
- The easy-fold system makes the iTravel the perfect travelling companion
- A walking stick holder is mounted as standard

OPERATION

- Very easy to use

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

ORTOPEDIA
EINE MARKE VON MEYRA
MEYRA QUALITY TESTED



Standard configuration

- Swing-up armrests
- Easy operation with Speed Preselection
- Battery charger
- Sensor
- Seat foldable
- Battery removable with additional charging socket

Allowed user weight

120 kg

Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply

Valid from 01.01.2026 – 31.12.2026.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

SCOOTER

Technical specifications

Length in mm	970
Width in mm	560
Seat depth in mm	400
Seat widths in mm	440
Height in mm	870
Back support height in mm	460
Height when folded in mm	300
Footplate to seat (lower leg length) in mm	440
Empty weight from	22 kg
User weight	120 kg
Driving speed	6 km/h
Castor wheels diameter	200
Drive wheels diameter	229
Ground clearance in mm	70
Max. obstacle height in mm	40
Turning radius in mm	970
Continuous motor performance	200 W
Peak motor performance	800 W
Batteries	10,4 Ah with lithium battery
Battery charger	2 A
Electronic	50 A
Empty weight with lithium batteries	23 kg
Permitted total weight in kg	150 kg

Please enter your quantities in the boxes

01 iTravel E-Mobile

Item	DESCRIPTION	Chassis colour	EUR €
9105414	MEYRA iTRAVEL – the travel wheelchair	White	3.637,00

02 Accessories

Item	DESCRIPTION	EUR €
10542837	Additional lithium battery	620,00