iTRAVEL CARBON W5905

1.074 Stock version

Allowed user weight

120 kg

MEYRA QUALITY TESTED



ORTOPEDIA

EINE MARKE VON MEYRA







Travel Wheelchair

ADAPTATION

- With its weight of 16 kg it is the lightweight among the electric wheelchairs
- The easy-fold system makes the iTravel the perfect travelling companion

OPERATION

- Very easy to use
- The armrest support folds back with
- Side panel with flip-up armrests for one-hand operation and easy transfer

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

ORTOPEDIA
EINE MARKE VON MEYRA
MEYRA QUALITY TESTED



Allowed user weight
120 kg

Hints
Our Terms and Conditions as stated at www.meyra.de/AGB shall apply

Valid from 01.01.2026 – 31.12.2026.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

SCOOTER

Technical specifications

Length in mm	923
Width in mm	570
Seat depth in mm	420
Seat widths in mm	450
Height in mm	930
Back support height in mm	460
Height when folded in mm	370
Footplate to seat (lower leg length) in mm	430
Empty weight from	16 kg
User weight	120 kg
Driving speed	6 km/h
Castor wheels diameter	178
Drive wheels diameter	215
Ground clearance in mm	70
Max. obstacle height in mm	40
Turning radius in mm	980
Continuous motor performance	*2x 250 W
Peak motor performance	1000 W
Batteries	12 Ah with lithium battery
Range	*8 km
Battery charger	2 A
Electronic	35 Ah
Empty weight with lithium batteries	17,8 kg
Permitted total weight in kg	140 kg

Please enter your quantities in the boxes

01 iTravel Carbon E-Mobile

	Item	DESCRIPTION	Chassis colour
	9107414	MEYRA iTRAVEL CARBON Lithium ion battery	Carbon