TA iQ IndoorWave

1.516

Allowed user weight

125 kg



MEYRA GROUP









with seat lift and tilt in space

CHARACTERISTICS

- Biomechanical backrest and footrest

FIELD OF USE

- For indoor and outdoor use
- The lowest seat height of 40 cm allows to drive under low objects and furniture
- Obstacle climbing up to 65 mm
- For independent mobility, for permanent and all-day use

THERAPEUTIC BENEFIT

- Particularly for users with more severe functional disorders such as high acquired or congenital spinal cord injury, cranial-brain trauma, advanced muscular diseases (e. g. muscular dystrophy), advanced MS sufferers, celebratory paresis, spina bifida

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please mail the completed order form to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

Order quantity:

MEYRA



Price basic model **15.608,00 €**
(without Back/Seating system)

Allowed user weight
125 kg

Hints
Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01.01.2026 – 31.12.2026.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

more product information at www.meyra.de

USER DATA

First name

Name

Street

ZIP CODE

City

Country

Birthday

Diagnosis 1

Diagnosis 2

Tel.

Mobil

E-mail

Order number

Health insurance fund / cost unit

Specialist trade Contact person/phone.

Supply specialist MEYRA

Approval info from week

ANATOMICAL DATA (IN MM)

1) Seat depth

2) Seat width

3) Chest width

4) Chest depth

5) Elbow height

6) Armpit height

7) Back height


8) Head height

9) Forearm length

10) LLL

Body height

Weight (kg)



Remarks

❑ I accept the Terms and Conditions of Meyra Distribution GmbH (valid as of January 2023)
https://www.meyra.com/gtc/

Date

City

Signature

TAS

Technical specifications

Driving speed	6 km/h
Gyroscope	not possible
Range	10 - 15 km
Max. User weight in kg	125 kg
Electric seat lift in mm	320
LED lighting	optional
Turning radius in mm	450
Batteries	2 x 38 Ah
Motor power	2 x 450 W
Electronic	R-Net
Length without leg supports in mm	695
Chassis width in mm	560
Seat height in mm	400
Seat widths in mm	380 / 400 / 450 / 500 / 550
Seat depth in mm	infinitely variable from 250 – 600
Max. obstacle height in mm	65
Ground clearance in mm	70
Fixing eyelets	optional
Electric seat angle adjustment in °	45°
Electric backrest adjustment	several options
Electric legrest	several options

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

01 Drive

	CODE	DESCRIPTION	RANGE	TAS	EUR €
❑	115	6 km/h-Variant	10 - 15 km		0,00

02 Wheels

		Pneumatic tyres		Puncture-proof tyres	
Castor	rim	Tire grey	Tire black	Tire grey	Tire black
Steering wheel 200 x 50	black	❑	❑	❑ 176,00	❑ 176,00
	silver	❑	❑	❑ 176,00	❑ 176,00

Drive wheels

	CODE	DESCRIPTION	TAS	EUR €
❑	4459	Drive wheel 260x85, grey PU tyres	95269	308,00
❑	4461	Drive wheel 260x85, black pneumatic tyres	95268	0,00

Type of drive

	CODE	DESCRIPTION	TAS	EUR €
❑	4444	Front wheel drive	95346	0,00

☐ Please select by ticking the relevant boxes (S) Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

03 Seat width / seat depth / seat frame

*ST (mm)	250	260	270	280	290	300	310	320	330	340	350	360	370	380	390	400	410	420	430	440	450	460	480	500	520	540	560	580	600
A											❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑
B	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑																		

	seat depth / seat frame		
Seat width (SW)	**Tension belt seat (SD mm)**	**Seat plate (SD mm)**	**Junior seat plate (SD mm)**
SW 250 mm			❑ B
SW 300 mm			❑ B
SW 350 mm			❑ B
SW 380 mm	❑ A	❑ A	
SW 400 mm	❑ A	❑ A	
SW 450 mm	❑ A	❑ A	
SW 500 mm	❑ A	❑ A	
SW 550 mm	❑ A	❑ A	

04 Back unit

	CODE	DESCRIPTION	TAS	EUR €
❑	401	Manual backrest 90°	80363	0,00
❑	4622	Manual backrest 90° Junior	81064	0,00
❑	26	Manual backrest, adjustable via gas spring	80577	544,00
❑	25	Electric backrest, 95°-135° ⚠ *Start, end position and adjustment range depend on seat depth and can vary by 15°*	97235	1.660,00
❑	4326	Electric backrest, 95°-150° ⚠ *Start, end position and adjustment range depend on seat depth and can vary by 15°*	80872/97234	3.267,00

⚠ *Back with physiological pivot*

TAS

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Varilite back system

	System width									
System	**250 mm**	**300 mm**	**350 mm**	**380 mm**	**400 mm**	**430 mm**	**450 mm**	**500 mm**	**550 mm**	**600 mm**
VL ICON-Back Deep	❑ 726,00		❑ 1.022,00	❑ 1.022,00	❑ 1.022,00	❑ 1.022,00	❑ 1.022,00	❑ 1.022,00		
VL ICON-Back Low			❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00		
VL ICON-Back Mid		❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00
VL ICON-Back Tall		❑ 1.068,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00		
without back system	❑									

05 Seat unit

	CODE	DESCRIPTION	TAS	EUR €
❑	27	Seat lift 320 mm		0,00
❑	118	Electrical seat angling (0 – 45°) *Compulsory coding*		3.012,00

without seat system

	CODE	DESCRIPTION	TAS	EUR €
❑	727	without seat system		0,00

☐ Please select by ticking the relevant boxes (S) Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

Back/seat options

	CODE	DESCRIPTION	TAS	EUR €
❑	913	Headrest with mounting	80337/80288	568,00
❑	4506	Headrest medium with mount	80337/80288	460,00
❑	4432	Headrest Junior with receptacle	80829/80288	460,00
❑	957	Torso pads (pair)	93569	168,00
❑	4336	Transport hook mounted on the back	80175	198,00
❑	4946	Lumbar support cushion small	93484	166,00
❑	4523	Lumbar support cushion medium	93484	166,00
❑	4524	Lumbar support cushion large	93484	166,00
❑	4525	Lumbar support cushion adjustable	93217	198,00
❑	292	Mesh fabric black		0,00
❑	4980	Incontinence material black		0,00

Thigh support

Position	rear	front
(A)	❑	❑

	CODE	FIXATION	DESCRIPTION	TAS	EUR €
❑	4333	(A)	Thigh supports (pair), 75 x 100 mm	80328/93262	874,00
❑	4334	(A)	Thigh supports (pair), 100 x 130 mm	80328/93263	874,00
❑	4335	(A)	Thigh supports (pair), 130 x 150 mm	80328/93264	874,00

Varilite Junior Seating System

		Cushion depth		
Varillte Junior	**CW**	**250 mm**	**300 mm**	**330 mm**
Single-chambered air cushion	250 mm	❑ 516,00		
	280 mm			❑ 516,00
	300 mm		❑ 516,00	

☐ Please select by ticking the relevant boxes (S) Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

Varilite Evolution seat system

		Cushion depth						
Varilite Evolution	**CW**	**350 mm**	**380 mm**	**400 mm**	**430 mm**	**450 mm**	**500 mm**	**550 mm**
Single-chambered air cushion	380 mm		❑ 1.068,00		❑ 1.068,00			
	400 mm			❑ 1.068,00		❑ 1.068,00		
	430 mm				❑ 1.068,00			
	450 mm			❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	500 mm					❑ 1.068,00	❑ 1.068,00	
	550 mm							

Varilite Meridian seat system

		Cushion depth						
Varilite Meridian	**CW**	**350 mm**	**380 mm**	**400 mm**	**430 mm**	**450 mm**	**500 mm**	**550 mm**
Dual-chambered air cushions	380 mm		❑ 1.068,00		❑ 1.068,00			
	400 mm	❑ 1.068,00		❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	430 mm				❑ 1.068,00			
	450 mm			❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	❑ 1.068,00
	500 mm			❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	550 mm					❑ 1.068,00	❑ 1.068,00	

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
(S) Special construction/individual solution o.r. - on request

Varilite Stratus seat system

		Cushion depth						
Varilite Stratus	**CW**	**350 mm**	**380 mm**	**400 mm**	**430 mm**	**450 mm**	**500 mm**	**550 mm**
Single-chambered air cushion	380 mm		❑ 624,00		❑ 624,00			
	400 mm			❑ 624,00		❑ 624,00		
	430 mm				❑ 624,00			
	450 mm			❑ 624,00		❑ 624,00	❑ 624,00	
	500 mm					❑ 624,00	❑ 624,00	

Varilite seating system Pro-Form NX

		Cushion depth							
Varilite Pro-Form NX	**CW**	**300 mm**	**350 mm**	**380 mm**	**400 mm**	**430 mm**	**450 mm**	**500 mm**	**550 mm**
Two chamber air cushion	300 mm		❑ 1.068,00						
	350 mm		❑ 1.068,00		❑ 1.068,00				
	380 mm			❑ 1.068,00		❑ 1.068,00			
	400 mm				❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	430 mm					❑ 1.068,00			
	450 mm				❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	500 mm				❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
(S) Special construction/individual solution o.r. - on request

06 Armrests swivel up

	with stabilizer (height 280-370mm)			without stabilizer (height 240-380mm)		
Guidance	**220 mm (jr.)**	**320 mm**	**430 mm**	**220 mm (jr.)**	**320 mm**	**430 mm**
Armrest left		❑	❑	❑ 1.162,00	❑ 1.162,00	❑ 1.162,00

	with stabilizer (height 280-370mm)			without stabilizer (height 240-380mm)		
Guidance	**220 mm (jr.)**	**320 mm**	**430 mm**	**220 mm (jr.)**	**320 mm**	**430 mm**
Armrest right		❑	❑	❑ 1.162,00	❑ 1.162,00	❑ 1.162,00

⚠ *Adjustable length and width position*

Removable armrests

	CODE	DESCRIPTION	TAS	EUR €
❑	106	Armrest at seat, height adjustable, removable	80347	744,00

Armrests options

	CODE	DESCRIPTION	TAS	EUR €
❑	4502	Adapter for lateral swivel armrests – right ⚠ *only in conjunction with CODE 14 swivel armrest*	80830	180,00
❑	4503	Adapter for lateral swivel armrests – left ⚠ *only in conjunction with CODE 14 swivel armrest*	80830	180,00
❑	4532	Protective edge for arm pads short	80651	73,00
❑	4533	Protective edge for arm pads long	80652	73,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

07 Legrests

Central legrests

Dimensions from seat plate

LLL (mm)	240	250	260	270	280	290	300	340	350	360	370	380	390	400	410	420	430
A								❑	❑	❑	❑	❑	❑	❑	❑	❑	❑
B	❑	❑	❑	❑	❑	❑	❑										

	CODE	LLL (mm)	DESCRIPTION	TYPE	TAS	EUR €
❑	5104	A	Multifunctional centra legrest incl. calf pad *with electric angle and lower leg length adjustment of 150mm*	Footplate	80373/80357	4.094,00
❑	5105	B	Multifunctional centra legrest incl. calf pad short *with electric angle and lower leg length adjustment of 80mm*	Footplate	80373/80571	4.094,00
❑	4510	A	Multifunctional centra legrest incl. calf pad *with electric angle and lower leg length adjustment of 150mm*	Single footplate	80818/81124	4.198,00
❑	4511	B	Multifunctional centra legrest incl. calf pad short *with electric angle and lower leg length adjustment of 80mm*	Single footplate	81126/80818	4.198,00
❑	4512	A	Multifunctional centra legrest incl. calf pad *with electric angle and lower leg length adjustment of 150mm*	electric footplate	81017/81128	6.772,00
❑	4513	B	Multifunctional centra legrest incl. calf pad short *with electric angle and lower leg length adjustment of 80mm*	electric footplate	81130/81017	6.772,00
❑	5101		Manual central legrest adjustable	Footplate	80373/95083	342,00
❑	5103		Manual central legrest adjustable	Single footplate	80818/95083	342,00
❑	4517		Manual central legrest adjustable	electric footplate	81133/95083	3.646,00
❑	4620		Manual legrest Central angle adjustable Junior	Footplate	80840/80841	342,00

Single legrest left

	CODE	DESCRIPTION	TAS	EUR €
❑	4450	Electric legrest	80629	1.942,00
❑	4451	Manual legrest	80701	262,00
❑	4508	Amputee leg support	80717	734,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Single legrest right

	CODE	DESCRIPTION	TAS	EUR €
❑	4330	Electric legrest	80628	1.942,00
❑	4331	Manual legrest	80700	262,00
❑	4507	Amputee leg support	80716	734,00

Leg support options

	CODE	DESCRIPTION	TAS	EUR €
❑	4329	Knee pads to support the central leg support ⚠ *only in conjunction with multifunctional legrest*	80690	358,00
❑	4498	Knee pad right side ⚠ *only in conjunction with CODE 4330 electric legrest*	80721	440,00
❑	4499	Knee pad left side ⚠ *only in conjunction with CODE 4450 electric legrest*	80730	440,00

08 Batteries and Chargers

	CODE	DESCRIPTION	TAS	EUR €
❑	509	Maintenance-free batteries 38 Ah	95542	0,00
❑	731	Charger 6 A	95388	0,00

09 Control modules and accessories

	CODE	DESCRIPTION	TAS	EUR €
❑	4462	R-Net LED control module	95174	0,00
❑	415	R-Net LED and light control module	95253	418,00
❑	416	R-Net control module CJSM2 and light	95252	0,00
❑	398	R-Net control module CJSM2 bluetooth	95524	636,00
❑	839	OMNI II Environment control for special controls ⚠ *Flashing function and horn function only selectable via menu*	80963/95541	2.834,00
❑	4789	Additional control module for accompanying person ⚠ *only in conjunction with CODE 4472, 4473 or 4467 bracket*	95213	1.118,00

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

Control module bracket

	CODE	DESCRIPTION	TAS	EUR €
❑	4327	Control module right, height adjustable, can be swiveled to the side	80407	266,00
❑	4328	Control module left, height adjustable, can be swiveled to the side	80408	266,00
❑	4504	Control module right, lowerable	80838	510,00
❑	4527	Control module left, lowerable	80838	510,00
❑	4472	Bracket for additional operating module ⚠ *only i.C. with Additional control Code 4789*	80517	308,00
❑	4473	Swivel bracket for additional operating module ⚠ *only i.C. with Additional control Code 4789*	80903	906,00

Control options

	CODE	DESCRIPTION	TAS	EUR €
❑	565	Chin control, electrically swing-away to the side	80692	5.514,00
❑	810	Chin control, mechanically swing-away to the side	80750	2.912,00
❑	4495	All-round Joystick ⚠ *only in conjunction with CODE 839 OMNI II*	80693	2.954,00
❑	4496	Switch Satellite Twister Black	80696	424,00
❑	4497	Multi-Joystick ⚠ *only in conjunction with CODE 839 OMNI II*	80704	5.808,00
❑	4500	Micro-Joystick ⚠ *only in conjunction with CODE 839 OMNI II*	80749	6.362,00
❑	4469	Heavy Duty Joystick ⚠ *only in conjunction with CODE 416 R-Net Joystick CJSM*	80839	3.356,00
❑	4505	All-round joystick, lightweight design ⚠ *only in conjunction with CODE 839 OMNI II*	81076	5.794,00
❑	4318	Multiswitch box (armrest) left or right	80678/97343	546,00
❑	4467	Hanger and emergency switch on the back for additional control module Code 4789 ⚠ *not in conjunction with transport hooks*	80533	2.970,00
❑	564	Rear-head control ⚠ *only in conjunction with CODE 839 OMNI II*	80951	6.782,00
❑	4337	Adapter for battery shutdown (transport)	95527	628,00

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

External button

	CODE	DESCRIPTION	TAS	EUR €
❑	4520	PIKO PB50, black	91003	326,00
❑	4521	PIKO PB50, black	91010	326,00
❑	4522	PIKO PB50, red	91010	326,00
❑	779	PIKO PB50, yellow	91010	340,00
❑	778	PIKO PB50, green	91010	340,00

Joystick options

	CODE	DESCRIPTION	TAS	EUR €
❑	4359	Softball 80 mm	95300	188,00
❑	4360	Golf ball white	80653	76,00
❑	4478	Golf ball black	80653	76,00
❑	4479	Golf ball pink	80653	76,00
❑	4480	Golf ball yellow	80653	76,00
❑	4481	Golf ball orange	80653	76,00
❑	4358	Tetra fork narrow 76 mm vertical	95301	228,00
❑	4593	Stielgriff 80 mm	95305	82,00
❑	4589	Mushroom handle 48mm	95154	70,00
❑	4588	Spherical handle 32mm	95299	70,00
❑	4590	Spherical handle 38mm	95155	70,00
❑	4594	Button handle 38mm	95306	70,00
❑	4586	Cone grip 30mm	95298	70,00
❑	4587	Joystick attachment standard	95153	0,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Mushroom, ball, knob and cone handles

	CODE	DESCRIPTION	TAS	EUR €
❑	4589	Mushroom handle 48mm		70,00
❑	4588	Ball handle 32 mm		70,00
❑	4590	Ball grip 38 mm		70,00
❑	4594	Button grip 38 mm		70,00
❑	4586	Cone grip 30mm		70,00

10 Lighting

	CODE	DESCRIPTION	TAS	EUR €
❑	932	LED lighting incl. turn signal, warning flashing light and brake light	97223	0,00
❑	703	Without lighting		0,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

11 Buckles for seat/back system

	CODE	DESCRIPTION	TAS	EUR €
		Breast shoulder strap		
❑	4401	XS	80765/93272	312,00
❑	4433	S	80765/93232	312,00
❑	4434	M	80765/93233	312,00
❑	4435	L	80765/93234	312,00
		Two-point belt		
❑	4436	XS	93221	188,00
❑	4437	S	93222	188,00
❑	4438	M	93223	188,00
❑	4439	L	93224	188,00
		Four- point belt		
❑	4440	XS	93225	220,00
❑	4441	S	93226	220,00
❑	4442	M	93227	220,00
❑	4443	L	93228	220,00

20 Accessories / Custom-made products

	CODE	DESCRIPTION	TAS	EUR €
❑	929	Therapeutic table transparent	80518/80163	1.338,00
❑	2721	Dahl Docking Adapter Plate ⚠ *Wheelchair adapter set for preparation of Dahl docking station*	95351	1.150,00
❑	721	Four-point fastening system for car transport	95350	0,00
❑	2237	Remote stop	80222	522,00
❑	953	Arm sleeve pocket	80994	612,00
❑	4518	Cladding black	95025	0,00
❑	4519	Cladding titanium	95027	106,00