TA iQ FWD STANDUP

1.520

Allowed user weight

140 kg

100 kg Function rising from lying (CODE 4464)

ISO-CRASH GETESTET



MEYRA GROUP









CHARACTERISTICS

- Fully suspensioned Powerbase (individual suspension concept for personal needs)
- Biomechanical backrest and footrest
- At 3 different speeds (6. 10. 13 km /h) available

FIELD OF USE

- For indoor and outdoor use
- The lowest seat height of 44 cm allows to drive under low objects and furnitures
- The seat and back are upholstered
- For independent mobility, for permanent and all-day use

THERAPEUTIC BENEFIT

- The individualised suspension concept, the fully-suspension chassis avoids unpleasant and painful vibrations and provides stability even on uneven terrain
- Especially for users with severe functional disorders, such as acquired or congenital tetraplegia, traumatic brain injury, advanced stage muscle diseases (e.g. muscular dystrophy), advanced stage MS sufferers, cerebral palsy, spina bifida
- The memory function enables recording of therapeutically motivated movement patterns as well as initial, end and preferred positions

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please mail the completed order form to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

Order quantity:

MEYRA



Price basic model **27.042,00 €**
(without Back/Seating system)

Standard configuration

- Front wheel drive
- Batteries
- Battery charger
- Seatlift
- Full suspension system
- R-NET Electronics
- Frame color black
- Gyroscope

Allowed user weight
140 kg
100 kg Function rising from lying (CODE 4464)

Hints
Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01.01.2026 – 31.12.2026.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

more product information at www.meyra.de

USER DATA

First name

Name

Street

ZIP CODE

City

Country

Birthday

Diagnosis 1

Diagnosis 2

Tel.

Mobil

E-mail

Order number

Health insurance fund / cost unit

Specialist trade Contact person/phone.

Supply specialist MEYRA

Approval info from week

ANATOMICAL DATA (IN MM)

1) Seat depth

2) Seat width

3) Chest width

4) Chest depth

5) Elbow height

6) Armpit height

7) Back height


8) Head height

9) Forearm length

10) LLL

Body height

Weight (kg)



Remarks

❑ I accept the Terms and Conditions of Meyra Distribution GmbH (valid as of January 2023)
https://www.meyra.com/gtc/

Date

City

Signature

TAS

Technical specifications

Driving speed	6 / 10 / 13 km/h
Gyroscope	Series
Range	30 - 40 km
Max. User weight in kg	140 kg (Rising from sitting) / 100 kg (Rising from lying)
Electric seat lift in mm	300
LED lighting	Series
Turning radius in mm	650
Batteries	2 x 80 Ah
Motor power	2 x 350 W
Electronic	R-Net
Length without leg supports in mm	880
Chassis width in mm	630
Seat height in mm	440
Seat widths in mm	380 / 400 / 450 / 500 / 550
Seat depth in mm	infinitely variable from 350 – 600
Max. obstacle height in mm	100
Ground clearance in mm	70
Fixing eyelets	Series
Electric seat angle adjustment in °	45°
Electric backrest adjustment	Series
Electric legrest	Series
Memory function	personalized
Combined function	Loungers / Seat stand

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

01 Drive

	CODE	DESCRIPTION	RANGE	TAS	EUR €
❑	115	6 km/h-Variant	40 km	99210/95380	0,00
❑	117	10 km/h-Variant	35 km	99210/95380	986,00
		⚠ *in Germany only in conjunction with CODE 4530 TüV certificate (category accessories)*			

02 Wheels

		Pneumatic tyres		Puncture-proof tyres	
Rear wheel	**rim**	**Tire grey**	**Tire black**	**Tire grey**	**Tire black**
Drive wheels 3.00-8	black	❑	❑	❑ 350,00	❑ 350,00
Drive wheels 3.00-8	silver	❑	❑	❑ 350,00	❑ 350,00

		Pneumatic tyres		Puncture-proof tyres	
Castor	**rim**	**Tire grey**	**Tire black**	**Tire grey**	**Tire black**
Steering wheel 200 x 50	black	❑	❑	❑ 176,00	❑ 176,00
	silver	❑	❑	❑ 176,00	❑ 176,00
Swivel wheel 210 x 65 reinforced forks	black	❑ 342,00	❑ 304,00	❑ 510,00	❑ 510,00

Wheels options

	CODE	DESCRIPTION	TAS	EUR €
❑	816	Plastic cover for swivel wheels 200 x 50 (pair)	97282	226,00
		⚠ *Not in conjunction with reinforced swivel wheels*		

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

03 Seat width/seat depth/seat frame

	*ST (mm)	380	390	400	410	420	430	440	450	460	480	500
A		❑	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑

Seat width (SW)	Seat depth (SD)	
SW 400 mm	❑	A
SW 450 mm	❑	A
SW 500 mm	❑	A

04 Varilite back system

System	System width 250 mm	300 mm	350 mm	380 mm	400 mm	430 mm	450 mm	500 mm	550 mm	600 mm
VL ICON-Back Deep			❑ 1.022,00	❑ 1.022,00	❑ 1.022,00	❑ 1.022,00	❑ 1.022,00	❑ 1.022,00		
VL ICON-Back Low			❑	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑		
VL ICON-Back Mid			❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00
VL ICON-Back Tall		❑ 1.068,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00	❑ 726,00		
without back system	❑									

TAS

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

05 Seat unit

	CODE	DESCRIPTION	TAS	EUR €
❑	4458	Electric stand-up function, including backrest, central legrest and armrests	99300	0,00
❑	4463	Biomechanical back length compensation		0,00
❑	4464	Combination function: lying in a standing position / sitting in a standing position ⚠ *max. user weight 100 kg*		0,00
❑	4465	Personalized memory function		0,00
❑	118	Electrical seat angling (0 – 45°) ⚠ *Compulsory coding*		2.868,00

without seat system

	CODE	DESCRIPTION	TAS	EUR €
❑	727	without seat system		0,00

☐ Please select by ticking the relevant boxes (S) Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

Back/seat options

	CODE	DESCRIPTION	TAS	EUR €
❑	913	Headrest with mounting	80337/80288	568,00
❑	4506	Headrest medium with mount	80337/80288	460,00
❑	957	Torso pads (pair)	93569	168,00
❑	4336	Transport hook mounted on the back	80175	198,00
❑	4946	Lumbar support cushion small	93484	166,00
❑	4523	Lumbar support cushion medium	93484	166,00
❑	4524	Lumbar support cushion large	93484	166,00
❑	4525	Lumbar support cushion adjustable	93217	198,00
❑	292	Mesh fabric black		0,00
❑	4980	Incontinence material black		0,00

Thigh support

Position	rear	front
A	❑	❑

	CODE	FIXATION	DESCRIPTION	TAS	EUR €
❑	4333	A	Thigh supports (pair), 75 x 100 mm	93262/80328	874,00
❑	4334	A	Thigh supports (pair), 100 x 130 mm	93263/80328	874,00
❑	4335	A	Thigh supports (pair), 130 x 150 mm	93264/80328	874,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Varilite Evolution seat system

		Cushion depth						
Varilite Evolution	**CW**	**350 mm**	**380 mm**	**400 mm**	**430 mm**	**450 mm**	**500 mm**	**550 mm**
Single-chambered air cushion	380 mm		❑ 1.068,00		❑ 1.068,00			
	400 mm			❑ 1.068,00		❑ 1.068,00		
	430 mm				❑ 1.068,00			
	450 mm			❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	500 mm					❑ 1.068,00	❑ 1.068,00	
	550 mm							

Varilite Meridian seat system

		Cushion depth						
Varilite Meridian	**CW**	**350 mm**	**380 mm**	**400 mm**	**430 mm**	**450 mm**	**500 mm**	**550 mm**
Dual-chambered air cushions	380 mm		❑ 1.068,00		❑ 1.068,00			
	400 mm	❑ 1.068,00		❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	430 mm				❑ 1.068,00			
	450 mm			❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	❑ 1.068,00
	500 mm			❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	550 mm					❑ 1.068,00	❑ 1.068,00	

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Varilite Stratus seat system

		Cushion depth						
Varilite Stratus	**CW**	**350 mm**	**380 mm**	**400 mm**	**430 mm**	**450 mm**	**500 mm**	**550 mm**
Single-chambered air cushion	380 mm		❑ 624,00		❑ 624,00			
	400 mm			❑ 624,00		❑ 624,00		
	430 mm				❑ 624,00			
	450 mm			❑ 624,00		❑ 624,00	❑ 624,00	
	500 mm					❑ 624,00	❑ 624,00	

Varilite seating system Pro-Form NX

		Cushion depth							
Varilite Pro-Form NX	**CW**	**300 mm**	**350 mm**	**380 mm**	**400 mm**	**430 mm**	**450 mm**	**500 mm**	**550 mm**
Two chamber air cushion	300 mm		❑ 1.068,00						
	350 mm		❑ 1.068,00		❑ 1.068,00				
	380 mm			❑ 1.068,00		❑ 1.068,00			
	400 mm				❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	430 mm					❑ 1.068,00			
	450 mm				❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	
	500 mm				❑ 1.068,00		❑ 1.068,00	❑ 1.068,00	

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

06 Armrest cushions

	CODE	DESCRIPTION	TAS	EUR €
❑	4322	Armrest cushion left, short 320 mm, length and width adjustable	80336	0,00
❑	4324	Armrest cushion left, length 430 mm, length and width adjustable	80321	0,00
❑	4323	Armrest cushion right, short 320 mm, length and width adjustable	80320	0,00
❑	4325	Armrest cushion right, length 430mm, length and width adjustable	80335	0,00

07 Legrests

Central legrests

Dimensions from seat plate

LLL (mm)	340	350	360	370	380	390	400	410	420	430
A	❑	❑	❑	❑	❑	❑	❑	❑	❑	❑

	CODE	LLL (mm)	DESCRIPTION	TYPE	TAS	EUR €
❑		A	Multifunctional centra legrest incl. calf pad	Footplate		0,00
			⚠ *with electric angle and lower leg length adjustment of 150mm*			

Leg support options

	CODE	DESCRIPTION	TAS	EUR €
❑	4531	Gel knee pads, set	81043	412,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

08 Batteries and Chargers

	CODE	DESCRIPTION	TAS	EUR €
❑	5158	Maintenance-free batteries 80 Ah	97332	0,00
❑	732	Charger 8 A	97236	0,00
❑	4316	Charger 10 A	97284	292,00

09 Control modules and accessories

	CODE	DESCRIPTION	TAS	EUR €
❑	4317	R-Net control module CJSM2	95501	398,00
❑	398	R-Net control module CJSM2 bluetooth	95524	636,00
❑	416	R-Net control module CJSM2 and light	95252	0,00
❑	839	OMNI II Environment control for special controls ⚠ *Flashing function and horn function only selectable via menu*	80963/95541	2.834,00
❑	4789	Additional control module for accompanying person ⚠ *only in conjunction with CODE 4472, 4473 or 4467 bracket*	95213	1.118,00

Control module bracket

	CODE	DESCRIPTION	TAS	EUR €
❑	4327	Control module right, height adjustable, can be swiveled to the side	80407	266,00
❑	4328	Control module left, height adjustable, can be swiveled to the side	80408	266,00
❑	4504	Control module right, lowerable	80838	510,00
❑	4527	Control module left, lowerable	80838	510,00
❑	4472	Bracket for additional operating module ⚠ *only i.C. with Additional control Code 4789*	80517	308,00
❑	4473	Swivel bracket for additional operating module ⚠ *only i.C. with Additional control Code 4789*	80903	906,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Control options

	CODE	DESCRIPTION	TAS	EUR €
❑	565	Chin control, electrically swing-away to the side	80692	5.514,00
❑	810	Chin control, mechanically swing-away to the side	80750	2.912,00
❑	4495	All-round Joystick ⚠ *only in conjunction with CODE 839 OMNI II*	80693	2.954,00
❑	4496	Switch Satellite Twister Black	80696	424,00
❑	4497	Multi-Joystick ⚠ *only in conjunction with CODE 839 OMNI II*	80704	5.808,00
❑	4500	Micro-Joystick ⚠ *only in conjunction with CODE 839 OMNI II*	80749	6.362,00
❑	4469	Heavy Duty Joystick ⚠ *only in conjunction with CODE 416 R-Net Joystick CJSM*	80839	3.356,00
❑	4505	All-round joystick, lightweight design ⚠ *only in conjunction with CODE 839 OMNI II*	81076	5.794,00
❑	4318	Multiswitch box (armrest) left or right	97343/80678	546,00
❑	564	Rear-head control ⚠ *only in conjunction with CODE 839 OMNI II*	80951	6.782,00
❑	4337	Adapter for battery shutdown (transport)	97285	628,00

External button

	CODE	DESCRIPTION	TAS	EUR €
❑	4520	PIKO PB50, black	91003	326,00
❑	4521	PIKO PB50, black	91010	326,00
❑	4522	PIKO PB50, red	91010	326,00
❑	779	PIKO PB50, yellow	91010	340,00
❑	778	PIKO PB50, green	91010	340,00

☐ Please select by ticking the relevant boxes (S) Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

Joystick options

	CODE	DESCRIPTION	TAS	EUR €
❑	4359	Softball 80 mm	95300	188,00
❑	4360	Golf ball white	80653	76,00
❑	4478	Golf ball black	80653	76,00
❑	4479	Golf ball pink	80653	76,00
❑	4480	Golf ball yellow	80653	76,00
❑	4481	Golf ball orange	80653	76,00
❑	4358	Tetra fork narrow 76 mm vertical	95301	228,00
❑	4593	Stielgriff 80 mm	95305	82,00
❑	4589	Mushroom handle 48mm	95154	70,00
❑	4588	Spherical handle 32mm	95299	70,00
❑	4590	Spherical handle 38mm	95155	70,00
❑	4594	Button handle 38mm	95306	70,00
❑	4586	Cone grip 30mm	95298	70,00
❑	4587	Joystick attachment standard	95153	0,00

10 Lighting

	CODE	DESCRIPTION	TAS	EUR €
❑	932	LED lighting incl. turn signal, warning flashing light and brake light	99116	0,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

11 Buckles for seat/back system

	CODE	DESCRIPTION	TAS	EUR €
		Breast shoulder strap		
❑	4401	XS	93272/80765	312,00
❑	4433	S	93232/80765	312,00
❑	4434	M	93233/80765	312,00
❑	4435	L	93234/80765	312,00
		Two-point belt		
❑	4436	XS	93221	188,00
❑	4437	S	93222	188,00
❑	4438	M	93223	188,00
❑	4439	L	93224	188,00
		Four- point belt		
❑	4440	XS	93225	220,00
❑	4441	S	93226	220,00
❑	4442	M	93227	220,00
❑	4443	L	93228	220,00

20 Accessories / Custom-made products

	CODE	DESCRIPTION	TAS	EUR €
❑	929	Therapeutic table transparent	80163	1.338,00
❑	2237	Remote stop	80222	522,00
❑	953	Armrest pocket	80994	612,00
❑	4530	TUV-certificate only i.c.w. code 117 ⚠ *only in conjunction with CODE 676 rear marking panel*		512,00
❑	2676	Warning sign sticker		284,00