iCHAIR MC2

1.611

Allowed user weight

140 kg with seat lift CODE 27 (incl. additional load)

160 kg (incl. additional load)



MEYRA GROUP



Electrically adjustable legrest available as an option



Comfortable operation



Active LED lighting



Castor wheels with aluminium rim

The intelligent and variable Allrounder

EFFICIENCY

- Effectiveness through simple maintenance and a high level of service friendliness
- Little need for repairs by implementing high value and durable parts
- Modular structure of the iChair-family enables compatibility of parts and parts groups
- Numerous adjustment options allow individual adaptation
- Simplified reinstallment through perfect user size limitations, as for example a standard user weight of 160 kg

THERAPEUTIC BENEFIT

- Ideal for retaining mobility at home, at the desk or in the office
- Effective pressure relief and relaxing positioning for long term sitting, through large electric adjustment sectors
- Different seat systems according to the individual demands

FIELD OF USE

- For indoor and outdoor use
- For permanent and all-day use for independent mobility, also with extreme functional limitations
- For working and active people, who need one wheelchair for different everyday situations
- Suitable for outdoor use through large serial tyres and optionally available 10 km/h-version
- Indoor spaces can easily be mastered due to the compact dimensions and high level of agility

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

Order quantity:



Allowed user weight

140 kg with seat lift CODE 27 (incl. additional load)
160 kg (incl. additional load)

Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01.01.2026 – 31.12.2026.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

USER DATA

First name

Name

Street

ZIP CODE

City

Country

Birthday

Diagnosis 1

Diagnosis 2

Tel.

Mobil

E-mail

Order number

Health insurance fund / cost unit

Specialist trade Contact person/phone.

Supply specialist MEYRA

Approval info from week

ANATOMICAL DATA (IN MM)

1) Seat depth

2) Seat width

3) Chest width

4) Chest depth

5) Elbow height

6) Armpit height

7) Back height


8) Head height

9) Forearm length

10) LLL

Body height

Weight (kg)



Remarks

❑ I accept the Terms and Conditions of Meyra Distribution GmbH (valid as of January 2023)
https://www.meyra.com/gtc/

Date City Signature

EPOWER

Technical specifications

Dimension tolerances +/- 15 mm, +/- 2°

1) *see category seat width*
2) *seat angle 4°*
3) *without legrests and without armrests*
4) *60 mm with wheel size 9"/12" and 70 mm with 10"/14", 120 mm incl. step climber*
5) *The adjustment range may deviate by max. +/-3° depending on SH and castor wheel size*
6) *with armrest CODE 106*
7) *up to height of seat pan*
8) *seating systems belt / ErgoSeat / ERGOSTAR / RECARO*
9) *9"/10"castor wheel*

Seat widths in mm	380 - 650 [1) 2)]
Seat depth in mm	400 - 560 / 400 - 560 / 500 / 480 - 530 [2) 8)]
Back support height in mm	450 - 500 / 530 - 570 / 640 / 640 [2) 8)]
Lower shank length in mm	280 - 430 [2) 7)]
Arm support height in mm	240 - 350 [7) 9)]
Front seat height in mm	440 - 530
Seat height with Code 27 seat lift	460 - 760 / 490 - 790
Length with footrest in mm	1130
Length without footrests in mm	830
Width ready to drive in mm from	590 - 630
Height without head support in mm	930 - 1130
Transport length in mm	830
Transport height in mm	570
Front wheel size in mm	230 x 70 (9") / 260 x 70 (10")
Rear wheel size in mm	320 x 60 (12,5") / 356 x 75 (14")
Max. obstacle height in mm	70 [4)]
Turning radius in mm	840
Permitted inclination/slopes	8,5° (15%)
Electrically adjustable backrest angle	80° - 140° (Code 25)
Mechanical backrest angle	80° - 120° (Code 401)

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

00 Chassis/Motorisation

	CODE	SPEED	RANGE	MOTOR POWER
❑	115	6 km/h	max. 40 km	2 x 220 W
❑	115/903	6 km/h ⚠ *Only in combination with R-Net*	max. 40 km	2 x 300 W (reinforced)
❑	117	10 km/h	max. 35 km	2 x 300 W
❑	117/903	10 km/h ⚠ *Only in combination with R-Net*	max. 35 km	2 x 350 W (reinforced)

ESP drive stabilisation module

	CODE	DESCRIPTION
❑	918	ESP drive stabilisation module (GYRO) ⚠ *Only in combination with R-Net LCD*

Chassis options

	CODE	DESCRIPTION
❑	892	Standard suspension for user weight up to 120 kg
❑	891	Reinforced suspension for user weight 120 - 160 kg

01 Seat width (SW)

	CODE	DESCRIPTION	SEAT SURFACE WIDTH
❑	38	SW continuously adjustable from 380 to 500 mm through the arm supports	380 mm
❑	43	SW continuously adjustable from 430 to 550 mm through the arm supports	430 mm
❑	48	SW continuously adjustable from 480 to 600 mm through the arm supports	480 mm
❑	53	SW continuously adjustable from 530 to 650 mm through the arm supports	530 mm

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

02 Front seat height (FSH) in combination with castor wheel (without seat cushion)

FSH measured at 4° seat inclination, upper edge of seat plate

		Front Seat Height (FSH) in mm					
Code	**Ø Rim**	**440**	**460**	**480**	**500**	**510**	**530**
163/146	9" aluminium puncture-proof pneumatic tyres	❑		❑		❑	
163/900	9" aluminium puncture-proof tyres	❑		❑		❑	
142/146	10" aluminium puncture-proof pneumatic tyres		❑		❑		❑
142/900	10" aluminium puncture-proof tyres		❑		❑		❑

⚠ *in conjunction with C. 27 and c. 4118, only VSH 510 mm for 9" (C. 163) or VSH 530 mm for 10" (C. 142) are possible.*

Steering wheel options

	CODE	DESCRIPTION
❑	816	Mud guard for Castor wheels

03 Drive wheels

		Puncture- resistant pneumatic tyres CODE 846		Puncture proof tyre CODE 904	
		(C-909)	(C-908)	(C-909)	(C-908)
Code	**Ø Rim**	**grey**	**black**	**grey**	**black**
815	12" aluminum rim	❑	❑	❑	❑
460	14" alloy rim	❑	❑	❑	❑

Rear wheel options

	CODE	DESCRIPTION
❑	824	Mudguard for drive wheels, Curved

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

05 Seat Tilt

	CODE	DESCRIPTION	Seat angle
❑	119	Mechanic seat tilt, in steps	0° to +10°
❑	118	Powered seat tilt continuously by 28 degrees ⚠ *Not possible with VR2 without adjustments*	-2° to +26°
❑	4118	Powered seat tilt stepless by 50 degrees ⚠ *Seat depth 56 cm CODE 356 is not possible in combination with ErgoSeat backrest CODE 949* *Only in conjunction with R-NET* *Seat height 460-760 mm with 9" castor wheels or 490-790 mm with 10" castor wheels (measured at 0° seat tilt)*	0° - 50°
❑	27	Scissor lift, 300 mm lift incl. seat tilt ⚠ *Only in conjunction with R-NET* *Seat height 460-760mm with 9" castor wheels or 490-790mm with 10" castor wheels (measured at 0° seat angle)*	0° to +30°

06 Seat depth (SD)

	(C-615)	(C-788)	(C-619)	(C-792)	(C-793)	(C-356)
	400 mm	430 mm	460 mm	490 mm	530 mm	560 mm
Seat depth (SD)	❑	❑	❑	❑	❑	❑

07 Seat system

	CODE	DESCRIPTION
❑	287/738	Seat plate without cushion
❑	287/561	Seat plate with cushion
❑	4276	Seat plate without cushion
❑	948	Ergo Seat cushion
❑	7206	Netti Uno Seat cushion
❑	7203	Netti Sit cushion

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

08 Back systems

	CODE	DESCRIPTION	BACK HEIGHT
❑	287/737	Standard Back (Canvas)	
❑	608	Round tube back posts, 22 mm, incl. continuous push bar ⚠ *prepared for external back systems*	
❑	735	Tension adjustable back strap without padding/back cover.	500 mm, convertible to 450 mm
❑	287/736	Tension adjustable back strap with padding / back cover.	500 mm, convertible to 450 mm
❑	949	ErgoSeat back,	590 mm, convertible to 630 mm
❑	2019	ErgoSeat Pro comfort back including overlay cushion	590 mm, convertible to 630 mm
❑	7201	Netti Uno Back rest	620 mm, convertible to 670 mm
❑	7209	Netti Smart Back rest	620 mm, convertible to 670 mm
❑	221/7207	Netti Stabil Back rest	620 mm
❑	638/7207	Netti Stabil Back rest	670 mm
❑	221/7210	Netti Super Stabil Back rest	620 mm
❑	638/7210	Netti Super Stabil Back rest	670 mm
❑	7205	Netti Kyphotic Back rest	620 mm, convertible to 670 mm

Optional seat units

Please consider the front height of the seat cushion: ERGOSTAR + 120 mm and RECARO + 150 mm

		Backrest	
CODE	Seat unit	mechanical	electric
961	ERGOSTAR anatomically formed seat system	❑	
579/567	RECARO S/F short seat with little side surface	❑	❑
580/578	RECARO M/W long seat with very long side bolsters	❑	❑
581/584	RECARO L/X long seat with pronounced lateral guidance	❑	❑

⚠ *ERGOSTAR seat CODE 961 is only possible with seat width 48 cm*

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

08 ErgoSeat-covers

	CODE	DESCRIPTION
❑	237	ErgoSeat cover, textile, black, with grey seams
❑	4980	ErgoSeat-cover, incontinenc, black

ErgoSeatPRO accessory back

	CODE	DESCRIPTION	Mass in mm (H X W Inside/Outside x D)
❑	2021	Thorax cushions, pair, large	310 x 125/70 x 140
❑	2032	Thorax cushions, pair, small	230 x 120/60 x 125

09 Recaro remuneration

	CODE	DESCRIPTION
❑	235	Recaro cover, textile, black

10 Back support options (Adjustments)

	CODE	DESCRIPTION
❑	814	Continuous push bar ⚠ *not in combination with RECARO or ERGOSTAR seating system*
❑	410	Single push handles incl. stabiliser bar
❑	2814	Through push bar for RECARO / ERGOSTAR seating system ⚠ *only in combination with RECARO or ERGOSTAR*
❑	913	Head support, synthetic leather black
❑	589	Shoulder Harness
❑	401	Mechanical back recline -10° - 30°/10° steps
❑	25	Powered back recline -10° to 50° stepless ⚠ *It is recommended to use the headrest, imitation leather black (CODE 913).*

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

11 Arm Supports and Arm Pads

1 = detachable 2 = height-adjustable 3 = width-adjustable 4 = depth-adjustable 5 = swivel-up 6 = angle-adjustable

	CODE	DESCRIPTION	FUNCTIONS
❑	106	Single post arm support with upholstered arm pad and black clothes guard	1,2,3,4
❑	4960	Single post arm support with upholstered arm support, satined clothes guard and slight depression in the padded arm support	1,2,3,4
❑	21	Single post arm support with upholstered arm support, satined clothes guard and long angled padded arm support	1,2,3,4
❑	24	Arm support with padded arm support, swivelling and angle adjustable Differing arm support height 190 - 230 mm ⚠ *Only i.c. with operating module bracket CODE 4859*	2,5,6
❑	14	Flip up armrests mounted on the back posts, adjustable in width, height and angle ⚠ *not in v.v.m. Code 954 and 957 Pads (pair)*	2,3,4,5,6
❑	769	Satined clothes guard without arm support ⚠ *only i.c. with arm support CODE 24*	
❑	4422	Armrest fixation in initial position ⚠ *only i.c.w. armrest CODE 14*	

Armrest hemiplegia

	CODE	DESCRIPTION
❑	4982	Hemiplegia forearm rest, swing-away a. height adjustable ⚠ *Not i.c.w. CODE 24 Side panel with padded armrest, swivel-up* *always opposite to the control module side* *For forearm length to the middle of the hand 36cm +/- 2cm*
❑	4983	Extension for hemiplegia forearm shell ⚠ *only i.c.w. CODE 4982 Hemi forearm support* *For forearm length to the middle of the hand 43cm +/- 2cm*

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

Hand Rest Hemiplegia

	CODE	DESCRIPTION
❑	4985	Hemi hand rest flat ⚠ *only in combination with CODE 4983 Hemi forearm tray*
❑	4987	Hemi hand rest with ergonomic handle on the left ⚠ *only in combination with CODE 4983 hemi forearm tray Control module side right CODE 60*
❑	4986	Hemi hand rest with ergonomic handle on the right ⚠ *only in combination with CODE 4983 hemi forearm tray Media module page left CODE 74*

12 Leg supports

	CODE	DESCRIPTION
❑	93	80° Legrests, detachable and swing-away
❑	92	Elevating Leg supports, (manual) , removable and swing-away ⚠ *not i.c.w.footrest CODE 54*
❑	86	Powered Elevating Legrests, removable and swing-away ⚠ *only i.c.w. operating module R-NET CODE 702* *not i.c.w. footrest CODE 54*

13 Footrests

	CODE	DESCRIPTION	depth	LLL
❑	54	One Piece footrest, 60 mm depth adjustable (in 4 steps), flip-up, height and angle adjustable ⚠ *Only i.c.w. legrest Code 93*	150 mm	280-430mm
❑	798	Two-piece footrest flip up, height and angle adjustable	150 mm	280-430mm
❑	105	Footboard adapter: abduction +/- 45 mm per side		
❑	822	Toe strap (pair)		
❑	823	Heel Loops (pair)		

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

14 Batteries and Chargers

	CODE	DESCRIPTION	RANGE
❑	5159	Maintenance free batteries approx. 43 Ah (5h), 50 Ah (20h)	25 km
❑	5158	Maintenance-free batteries approx.66 Ah (5h), 80 Ah (20h)	40 km
❑	732	Charger 6 A for 50 Ah(20h)	
❑	733	Charger 12A recommended for 80 Ah (20 h)	

15 Lighting

	CODE	DESCRIPTION
❑	932	Active LED-lighting front and rear acc. to StVZO

16 Joystick modules and accessories

		CONTROL			
Code	Control modules	none	2 / 682	6 / 678	8 / 780
408	VR2	❑	❑		
702/415	R-NET LED			❑	
702/416	R-NET LCD			❑	
702/398	R-NET LCD IR/BT			❑	
702/839	R-NET OMNI II			❑	

⚠ *Code 839 OMNI II Flashing function and horn function only selectable via menu*

Joystick Location

	CODE	DESCRIPTION	CONTROLS
❑	74	Assembly right	
❑	60	Assembly left	

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

17 Joystick Mount

	CODE	DESCRIPTION
❑	4859	Joystick bracket, length adjustable
❑	4566	Joystick bracket, length and height adjustable
❑	4860	Joystick bracket, length and height adjustable, swivelling to the side (e.g. on a table)
❑	4874	Protection bar for Joystick VR2 and R-NET
❑	851	Joystick bracket for accompanying person (operating module can be moved from the front to the rear) ⚠ *Not compatible with single push handles including stabilising bar CODE 410*

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

18 Special operating modules

	CODE	DESCRIPTION
❑	4660	Table plate with recess (centred) for operating module, swivelling to the right
❑	4661	Table plate with recess (centred) for operating module, swivelling to the left
❑	4789	Dual Control Joystick with priority switch for accompanying person, incl. assembly on the frame tube
❑	596	External on/off switch, diameter 63 mm, yellow ⚠ *Only i.c.w. joystick module R-NET with LED-Display CODE 702/416 or OMNI II Surround control CODE 839*
❑	778	External on/off switch, diameter 30 mm, green ⚠ *Only i.c.w. joystick module R-NET with LED-Display CODE 702/416 or OMNI II Surround control CODE 839*
❑	779	External on/off switch, diameter 30 mm, yellow ⚠ *Only i.c.w. joystick module R-NET with LED-Display CODE 702/416 or OMNI II Surround control CODE 839*

Zero displacement joystick

	CODE	DESCRIPTION
❑	4940	Neck holder bracket ⚠ *Only in conjunction with CODE 4611 or 4612*
❑	4612	MicroGuide way joystick ⚠ *only i. c. with R-NET OMNI II CODE 702/839* *temporarily not available*
❑	4601	MicroGuide way joystick including bracket and hand pad ⚠ *only in combination with R-NET OMNI II CODE 702/839* *temporarily not available*

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Joystick attachments

	CODE	DESCRIPTION
❑	553	Tetra-fork vertical
❑	556	Tetra-fork horizontal
❑	734	Tetra-fork horizontal, adjustable to hand width
❑	4589	Mushroom grip
❑	4590	Ball grip
❑	4591	T-grip
❑	4593	Stem grip
❑	4594	Button grip
❑	4595	Softball
❑	4596	Golf ball d=43 mm

19 Frame colours / accent colours

	CODE	FRAME	AKZENTE
❑	207	Matt black	
❑	212	Matt black	Polar white
❑	186	Matt black	Magic Blue
❑	227	Matt black	Red metallic
❑	180	Matt black	Lemon reflex
❑	214	Matt black	Sun yellow

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

20 Options around the wheelchair

	CODE	DESCRIPTION
❑	906	Rear-view mirror right
❑	927	Rear-view mirror left
❑	994	Baggage rack with metal rods
❑	970	Crutch holder
❑	4795	Anti-theft function, key ⚠ *only i.c. with R-Net*
❑	833	Safety belt (metal buckle)
❑	721	Lap belt - buckle
❑	4646	Therapy table swivelling sideways
❑	590	Abduction wedge ⚠ *Only i.c.w. ErgoSeat CODE 948 or seat plate CODE 4276*
❑	862	Kerb climber
❑	260	Memory-function for head support
❑	957	Lateral / Torso supports, Rigid(pair) ⚠ *not in conjunction with side panel CODE 14*
❑	954	Lateral / Torso supports, swing away (pair) ⚠ *not in conjunction with side panel CODE 24 not in conjunction with side panel CODE 14*
❑	781	USB charging socket
❑	782	Receptacle for mobile phone bracket
❑	2721	Adapter set for DAHL docking station
❑	2676	Warning sign sticker
❑	676	Rear marking plate

21 Service

	CODE	Item	DESCRIPTION
❑		55704	Flat fee for customised parameterisation (not eligible for discounts)

⚠ *The on-the-spot service covers only simple adjustment of parameter sets, such as acceleration, speed etc.*