iCHAIR ORBIT

1.618

Allowed user weight

120 kg

160 kg (with 350 watt motors CODE 525)



MEYRA GROUP



Optimal maneuverability and high stability for precise maneuvering



Compact dimensions realized by the smallest overall width and seat height



Maximum flexibility realized by single or double support wheel



High max. user weight of 160 kg

The flexible agile one

EFFICIENCY

- Effectiveness through simple maintenance and a high level of service friendliness
- Little need for repairs by implementing high value and durable parts
- Modular structure of the iChair-family enables compatibility of parts and parts groups
- Numerous adjustment possibilities enable individual adaptation
- Simplified reinstallment through perfect user size limitations, as for example a standard user weight of 160 kg

THERAPEUTIC BENEFIT

- The central leg support makes a 90° knee angle possible
- Ideal for retaining mobility at home, at the desk or in the office
- Effective pressure relief and relaxing positioning for long term sitting, through large electric adjustment sectors

FIELD OF USE

- For indoor and outdoor use, spring-mounted drive unit standard
- Especially good for indoor use: Rotating around its own axle supports in tight manoeuvres
- Precise manoeuvres through rigidly mounted steering castors
- The optional side-mounted support wheels allow easy and safe driving on single ramps
- For permanent and all-day use for independent mobility, also with extreme functional limitations
- For working and active people, who need one wheelchair for different everyday situations

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

Order quantity:



Price basic model ***7.373,00 €***

Allowed user weight
120 kg
160 kg (with 350 watt motors CODE 525)

Hints
Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01.01.2026 – 31.12.2026.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

EPOWER

USER DATA

First name

Name

Street

ZIP CODE

City

Country

Birthday

Diagnosis 1

Diagnosis 2

Tel.

Mobil

E-mail

Order number

Health insurance fund / cost unit

Specialist trade Contact person/phone.

Supply specialist MEYRA

Approval info from week

ANATOMICAL DATA (IN MM)

1) Seat depth

2) Seat width

3) Chest width

4) Chest depth

5) Elbow height

6) Armpit height

7) Back height


8) Head height

9) Forearm length

10) LLL

Body height

Weight (kg)



Remarks

❑ I accept the Terms and Conditions of Meyra Distribution GmbH (valid as of January 2023)
https://www.meyra.com/gtc/

Date

City

Signature

EPOWER

Technical specifications

Dimension tolerances +/- 15 mm, +/- 2°
Data subject to design changes
1) *see category seat width*
2) *to seat pan without cushion*
3) *with adjustable back*
4) *seating Systems: Standard / ErgoSeat / Recaro*
5) *with 350W motors*
6) *measured at 0° seat angle, 12"/14"*
7) *Single/central legrests*
8) *CODE 38/43/48/53, with 350W + 10 mm total width*

Seat widths in mm	380 - 650 [1]
Seat depth in mm	400-560/400-560/480-530 [1]
Front seat height in mm	430, 465, 500/440, 475, 510 [6]
Back support height in mm	450 - 500/530 - 580/640 [4]
Lower shank length in mm	280 - 450 [2]
Arm support height in mm	240 - 350 [2]
Length with central footrests in mm	1180 - 1250/1060 - 1150 [7]
Length without footrests in mm	930
Width 12 inches ready to drive in mm	600/600/630/680 [8]
Width 14 inches ready to drive in mm	640/640/640/680 [8]
Height without head support in mm	980 - 1060
Transport height min. in mm	630
Max. obstacle height in mm	60
Turning radius with central legrest in mm	700
Turn radius with single leg support in mm	800
Ground clearance	70
Transport empty weight basic configuration min.	120 kg - 165 kg
Permitted overall weight incl. options	320 kg [5]
Max. additional loading	10 kg
Driving speed	6 km/h / 10 km/h
Motor power	2 x 180/2 x 300
Range with battery 50 Ah(20h)	25 km - 35 km
Range with battery 73 Ah (20h)	35 km - 40 km
Permitted inclination/slopes	8,5° (15%)
Form of steering	indirect
Seat camber electric	0° to +30°/0° to +50°
Seat camber mechanical	0°, 4°, 7°, 10°
Stand-up aid (0° to -25°) with electric tilt, lift 200 mm	Seat tilt/stroke -25° to 45
Back support angle electric	80° - 140° (Code 25)

Technical data *** continued

Back support angle mechanical	80° - 120° (Code 401)
Lift in mm	300
Drive	Midwheel

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

00 Chassis / engine / frame options

	CODE	SPEED	RANGE	MOTOR POWER	MAX. Weight.	EUR €
❑	115	6 km/h	max. 40 km	2 x 180 W	120 kg	0,00
		⚠ *only in combination with operating module VR2 CODE 408*				
❑	115/525	6 km/h	max. 40 km	2 x 350 W	160 kg	1.117,00
❑	117/525	10 km/h	max. 35 km	2 x 350 W	160 kg	1.585,00
		⚠ *in Germany only i.c. with CODE 4772 TÜV certificate (accessories section)* *in Germany only i.c with CODE 676 rear marking plate*				
❑	918	ESP drive stabilisation module (GYRO)				1.665,00
		⚠ *ESP driving stabilization module, only in conjunction with R-Net / LCD color display or i.V.m. OMNI II (CODE 839)*				

Frame design and support wheel options

	CODE	DESCRIPTION	EUR €
❑	352	Rear swing arm with central wheel suspension	0,00
❑	690	Frame with rear swing arm, LH/RH	176,00
❑	2001	Frame design wide	0,00

01 Seat width (SW)

	CODE	DESCRIPTION	SEAT SURFACE WIDTH	EUR €
❑	38	SW continuously adjustable from 380 to 440 mm through the arm supports	380 mm	0,00
❑	43	SW continuously adjustable from 430 to 520 mm through the arm supports	430 mm	0,00
❑	48	SW continuously adjustable from 480 to 600 mm through the arm supports	480 mm	0,00
❑	53	SW continuously adjustable from 530 to 650 mm through the arm supports	530 mm	0,00

⚠ *I.c.w. ErgoSeat Pro only SW 430 and 480 mm possible*
With RECARO the real seat width is approx. 500 mm
Overview of RECARO seats on request

02 Castor

Code	Ø	rim	Tyres	
157/146	8"	Aluminium rim	Punctured tires	❑
157/900	8"	Aluminium rim	Puncture-proof tyres	❑ 189,00

EPOWER

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Damped steering wheels

	CODE	DESCRIPTION	EUR €
❑	2145	8" cushioned castor forks (recommended for users weighing up to 100 kg)	944,00
❑	2146	8" cushioned castor forks (recommended for users weighing up to 160 kg)	944,00

⚠ *Only compatible with 14" drive wheel with aluminium rim CODE 460. Not compatible with Dahl docking station adapter set CODE 2721. Seat height*

03 Drive wheels in combination with front seat height (FSH)

FSH (mm)	430	440	465	475	500	510
A		❑		❑		❑
B	❑		❑		❑	

		Color tires (front seat height)		
CODE	Drive wheel / tires	black (abrasion resistant)	grey	black (cleat profile)
460/846	14" drive puncture-resistant pneumatic tyres with aluminium rims	❑ A 250,00	❑ A 250,00	❑ A 280,00
460/904	14" drive puncture-proof tyres wheels with aluminium rims	❑ A 447,00	❑ A 447,00	
815/846	12" aluminium rims, puncture-resistant pneumatic tyres	❑ B	❑ B	
815/904	12" aluminium rims, puncture-proof tyres	❑ B 197,00	❑ B 197,00	

⚠ *VSH measured at 0° seat tilt, up to seat plate height*

Rear wheel options

	CODE	DESCRIPTION	EUR €
❑	824	Mudguard for drive wheels, Curved	126,00

EPOWER

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

05 Seat Tilt

	CODE	DESCRIPTION	Seat angle	LIFT	EUR €
❑	119	Mechanical seat tilting	0°/4°		0,00
❑	118	Electric seat tilting ⚠ *only i.c. with adjustment module VR2 CODE 682 or R-Net CODE 702*	0° to 30°		1.424,00
❑	4118	Electric seat tilting ⚠ *The seat depth 56 cm CODE 356 is not possible in combination with ErgoSeat backrest CODE 949* *Only in conjunction with adjustment module VR2 CODE 682 or R-Net CODE 702*	0° to 50°		1.907,00
❑	27	Scissor-lift 300 mm incl. tilt module ⚠ *only i. c. with R-Net* *User weight up to 140 kg*	0° to 30°	300 mm	5.408,00
❑	4937	Stand-up aid (0° to -25°) with electric tilt for easy transfer from wheelchair ⚠ *only i. c. with CODE 4260 External adjustment module* *only i. c. with CODE 900 PU tyres* *recommended with CODE 5025 electric backrest adjustment with memory function* *Seat height can vary by up to 20 mm*	0° bis 45°	200 mm	5.771,00

06 Seat depth (SD)

	(C-615)	(C-788)	(C-619)	(C-792)	(C-793)	(C-356)
	400 mm	430 mm	460 mm	490 mm	530 mm	560 mm
Seat depth (SD)	❑	❑	❑	❑	❑	❑ 44,00

07 Seat system

	CODE	DESCRIPTION	EUR €
❑	561	Seat cushion	247,00
❑	4276	Seat plate without cushion	104,00
❑	948	Ergo Seat cushion	384,00
❑	7206	Netti Uno Seat cushion	344,00
❑	7203	Netti Sit cushion	564,00

⚠ *Netti and Ride not i.c. with SD 560 mm*

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

08 Back systems

	CODE	DESCRIPTION	BACK HEIGHT	EUR €
❑	608	Round tube back posts, 22 mm, incl. continuous push bar ⚠ *Only in combination with mechanical backrest adjustment CODE 401 or electric backrest adjustment CODE 25* *Prepared for external backrest systems*		0,00
❑	287/735	Tension adjustable back strap without padding/back cover.	500 mm, convertible to 450 mm	81,00
❑	287/736	Tension adjustable back strap with padding / back cover.	500 mm, convertible to 450 mm	132,00
❑	949	ErgoSeat back,	590 mm, convertible to 630 mm	374,00
❑	2019	ErgoSeat Pro comfort back including overlay cushion	590 mm, convertible to 630 mm	525,00
❑	7201	Netti Uno Back rest	620 mm, convertible to 670 mm	477,00
❑	7209	Netti Smart Back rest	620 mm, convertible to 670 mm	477,00
❑	221/7207	Netti Stabil Back rest	620 mm	729,00
❑	638/7207	Netti Stabil Back rest	670 mm	729,00
❑	221/7210	Netti Super Stabil Back rest	620 mm	825,00
❑	638/7210	Netti Super Stabil Back rest	670 mm	825,00
❑	7205	Netti Kyphotic Back rest	620 mm, convertible to 670 mm	963,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

09 Optional seat units

Please consider the front height of the seat cushion: ERGOSTAR + 120 mm and RECARO + 150 mm

		Backrest	
CODE	Seat unit	mechanical	electric
961	ERGOSTAR anatomically formed seat system	❑	
579/567	RECARO S/F short seat with little side surface	❑ 2.764,00	❑ 3.698,00
580/578	RECARO M/W long seat with very long side bolsters	❑ 2.522,00	❑ 4.299,00
581/584	RECARO L/X long seat with pronounced lateral guidance	❑ 2.921,00	❑ 4.299,00

ERGOSTAR seat CODE 961 is only possible with seat width 48 cm
Electric backrest adjustment only in conjunction with R-Net control
Longer delivery time for RECARO seats

10 ErgoSeat-covers

	CODE	DESCRIPTION	EUR €
❑	237	ErgoSeat cover, textile, black, with grey seams	0,00
❑	4980	ErgoSeat-cover, incontinenc, black	0,00

Options of the ErgoSeat Pro comfort back

	CODE	DESCRIPTION	Mass in mm	EUR €
❑	2021	Thorax cushions, pair, large	310 x 125/70 x 140	311,00
❑	2032	Thorax cushions, pair, small	230 x 120/60 x 125	294,00

11 Recaro remuneration

	CODE	DESCRIPTION	EUR €
❑	235	Recaro cover, textile, black	0,00

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

12 Back support options (Adjustments)

	CODE	DESCRIPTION	EUR €
❑	814	Continuous push bar	0,00
❑	410	Single push handles incl. stabiliser bar	145,00
❑	2814	Continuous push bar for RECARO seat systems ⚠ *only in combination with RECARO*	229,00
❑	913	Head support, synthetic leather black	269,00
❑	5025	Electric back adjustment with memory function only for stand-up aid code 4937 ⚠ *It is recommended to use the headrest, imitation leather black (CODE 913).*	1.579,00
❑	401	Mechanical backrest recline	0,00
❑	25	Powered backrest recline ⚠ *Only in conjunction with VR2 CODE 682 or R-Net CODE 702 adjustment module. It is recommended to use the headrest, imitation lea*	1.462,00

13 Arm Supports and Arm Pads

1 = detachable 2 = height adjustable 3 = width adjustable 4 = depth adjustable

	CODE	DESCRIPTION	FUNCTIONS	EUR €
❑	106	Single post arm support with upholstered arm pad and black clothes guard	1,2,3,4	0,00
❑	4960	Side part with padded armrest, black clothes guard with colour inlay and slight depression in the armrest	1,2,3,4	74,00
❑	21	Single post arm support with upholstered arm support, satined clothes guard and long angled padded arm support	1,2,3,4	340,00
❑	24	Side part with upholstered armrest, swiveling up and angle-adjustable ⚠ *only in combination with CODE 4859 operating module holder*	2	195,00
❑	14	Flip up armrests mounted on the back posts, adjustable in width, height and angle		1.063,00
❑	4422	Armrest fixation in initial position ⚠ *only i.c.w. armrest CODE 14*		59,00

Armrest hemiplegia

	CODE	DESCRIPTION	EUR €
❑	4982	Hemiplegia forearm rest, swing-away a. height adjustable ⚠ *Not i.c.w. CODE 24 Side panel with padded armrest, swivel-up* *always opposite to the control module side* *For forearm length to the middle of the hand 36cm +/- 2cm*	386,00
❑	4983	Extension for hemiplegia forearm shell ⚠ *only i.c.w. CODE 4982 Hemi forearm support* *For forearm length to the middle of the hand 43cm +/- 2cm*	53,00

EPOWER

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Hand Rest Hemiplegia

	CODE	DESCRIPTION	EUR €
❑	4985	Hemi hand rest flat ⚠ *only in combination with CODE 4983 Hemi forearm tray*	58,00
❑	4987	Hemi hand rest with ergonomic handle on the left ⚠ *only in combination with CODE 4983 hemi forearm tray Control module side right CODE 60*	58,00
❑	4986	Hemi hand rest with ergonomic handle on the right ⚠ *only in combination with CODE 4983 hemi forearm tray Media module page left CODE 74*	58,00

14 Leg supports

	CODE	DESCRIPTION	EUR €
❑	93	80° Legrests, detachable and swing-away	0,00
❑	92	Elevating Leg supports, (manual) , removable and swing-away ⚠ *not i.c.w.footrest CODE 54*	252,00
❑	86	Powered Elevating Legrests, removable and swing-away ⚠ *only i.c.w. operating module R-NET CODE 702* *not i.c.w. footrest CODE 54*	1.829,00

Central legrest

	CODE	DESCRIPTION	LS-ANGLE	LLL	EUR €
❑	5100	Central legrest, mechanically angle adjustable	85°-95°	280-400 mm	470,00
❑	5101	Central legrest, mechanically angle adjustable	85°-95°	330-450 mm	470,00
❑	4935	Central legrest electrically angle adjustable ⚠ *For VSH up to 440 mm* *Only in conjunction with CODE 702/678 R-Net with adjustment module* *Not i.c. with RECARO* *Only i.c. with CODE 945 Calf padding*	90°-180°	200-400mm	2.381,00
❑	5104	Central legrest electrically angle adjustable ⚠ *For FSH from 475 mm* *Only i.c.w. CODE 702/678 R-Net with adjustment module* *Not i.c.w RECARO* *Only i.c.w. CODE 945 Calf pad*	90°-180°	200-450mm	2.381,00
❑	4949	Central legrest, electrically angle/length adjustable, electrically lowerable to the floor, ergonomic length compensation up to LS angle 150°. ⚠ *Only i.c.w. CODE 702/678 R-Net with adjustment module* *Not i.c.w. RECARO* *Only i.c.w. CODE 945 calf pad*	90°-180°	350-450mm	2.915,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Footplate options for central legrest

	CODE	DESCRIPTION	EUR €
❑	5109	Standard footplate	0,00
❑	5108	large footplate ⚠ *only i.c. with code 2001 wide frame version*	0,00
❑	5107	Single footplates, separately adjustable ⚠ *only i.c. with code 2001 wide frame* *not i.c. with code 4949*	0,00

15 Footrests

	CODE	DESCRIPTION	depth	LLL	EUR €
❑	54	One Piece footrest, 60 mm depth adjustable (in 4 steps), flip-up, height and angle adjustable ⚠ *Only i.c.w. legrest Code 93*	150 mm	280-430mm	395,00
❑	798	Two-piece footrest flip up, height and angle adjustable	150 mm	280-430mm	0,00

Footrest options

	CODE	DESCRIPTION	EUR €
❑	105	Footboard adapter: abduction +/- 25 mm each side ⚠ *Not i.c.w. with central legrest*	120,00
❑	822	Toe strap (pair)	59,00
❑	823	Heel Loops (pair)	105,00
❑	2945	Calf cushion i.c. with mechanical central footrest	380,00

16 Batteries and Chargers

	CODE	DESCRIPTION	EUR €
❑	5159	Maintenance-free batteries	0,00
❑	5158	Maintenance-free batteries	290,00
❑	732	*Ladegerät 6A für 50 Ah	0,00
❑	733	Charger 12 A for 80 Ah	305,00

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

17 Lighting

	CODE	DESCRIPTION	EUR €
❑	932	Active LED-lighting front and rear acc. to StVZO	0,00

18 Joystick modules and accessories

		CONTROL			
Code	Control modules	none	2 / 682	6 / 678	8 / 780
408	VR2	❑	❑ 368,00		
702/415	R-NET LED			❑ 974,00	❑ 2.505,00
702/416	R-NET LCD			❑ 1.342,00	❑ 2.873,00
702/398	R-NET LCD IR/BT			❑ 2.150,00	❑ 3.681,00
702/839	R-NET OMNI II			❑ 3.320,00	❑ 4.851,00

2-fold control only for electr. seat tilt and electr. backrest adjustment
Code 839 OMNI II flashing function and horn function only selectable via menu

Joystick Location

	CODE	DESCRIPTION	EUR €
❑	74	Assembly right	0,00
❑	60	Assembly left	0,00

19 Joystick Mount

	CODE	DESCRIPTION	EUR €
❑	4859	Joystick bracket, length adjustable	0,00
❑	4860	Joystick bracket, length and height adjustable, swivelling to the side (e.g. on a table)	281,00
❑	4874	Protection bar for Joystick VR2 and R-NET	65,00

☐ Please select by ticking the relevant boxes (S) Special construction/individual solution o.r. - on request

Standard version: If no other selection is made, the standard version will be delivered. 351,00

⚠ *Only i.c.w. continuous pushing bar CODE 814*

662,00

⚠ *Not i.c.w. push handles CODE 410*

20 Special operating modules

Special operating modules and special controls often have to be adapted and can therefore generally not represent a finished product.
Placement and programming must be carried out with the involvement of the user or worked out by a technician on site.
Other special controls are possible on request.

	CODE	DESCRIPTION	EUR €
❑	4660	Table plate with recess (centred) for operating module, swivelling to the right ⚠ *Not i.c.w. CODE 24 side panel* *Not i.c.w. CODE 702/839 R-NET OMNI II*	710,00
❑	4661	Table plate with recess (centred) for operating module, swivelling to the left ⚠ *Not i.c.w. CODE 24 side panel* *Not i.c.w. CODE 702/839 R-NET OMNI II*	710,00
❑	4646	Therapy table swivelling sideways	649,00

External on/off switch

for users who can hardly operate the keys of the operating module or as an emergency/off function for accompanying person

	CODE	DESCRIPTION	EUR €
❑	596	External on/off switch, diameter 63 mm, yellow	313,00
❑	778	External on/off switch, diameter 30 mm, green	281,00
❑	779	External on/off switch, diameter 30 mm, yellow	265,00

⚠ *Only i.c.w. joystick module R-NET with LED-Display CODE 702/416 or OMNI II Surround control CODE 839*

Others

	CODE	DESCRIPTION	EUR €
❑	4612	MicroGuide way joystick ⚠ *only in combination with R-NET OMNI II CODE 702/839* *temporarily not available*	5.267,00
❑	4601	MicroGuide way joystick including bracket and hand pad ⚠ *only in combination with R-NET OMNI II CODE 702/839* *temporarily not available*	5.737,00

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

Joystick attachments

	CODE	DESCRIPTION	EUR €
❑	553	Tetra-fork vertical	204,00
❑	556	Tetra-fork horizontal	204,00
❑	734	Tetra-fork horizontal, adjustable to hand width	134,00
❑	4589	Mushroom grip	64,00
❑	4590	Ball grip	64,00
❑	4591	T-grip	126,00
❑	4593	Stem grip	69,00
❑	4594	Button grip	64,00
❑	4595	Softball	64,00
❑	4596	Golf ball, small (d=43 mm)	99,00

21 software

	CODE	DESCRIPTION	EUR €
❑	985	Series programming (standard)	0,00
❑	998	Special programming (deviation from standard)	0,00

22 Frame Colours

	CODE	DESCRIPTION	EUR €
❑	207	matt black structure (without rim inserts)	0,00
❑	212	polar white	120,00
❑	186	magic blue	120,00
❑	227	Redmetallic	120,00
❑	180	lemon reflex	120,00
❑	214	sun yellow	120,00

☐ Please select by ticking the relevant boxes
S Special construction/individual solution o.r. - on request
Standard version: If no other selection is made, the standard version will be delivered.

23 Options around the wheelchair

	CODE	DESCRIPTION	EUR €
❑	906	Rear-view mirror right	95,00
❑	927	Rear-view mirror left	95,00
❑	970	Crutch holder	95,00
❑	833	Safety belt (metal buckle)	91,00
❑	4646	Therapy table swivelling sideways	649,00
❑	590	Abduction wedge ⚠ *Not i.c.w. Central leg support CODE 5100/5101/5104/4949/4935*	275,00
❑	260	Memory-function for head support	56,00
❑	957	Lateral / Torso supports, Rigid(pair)	231,00
❑	954	Lateral / Torso supports, swing away (pair) ⚠ *not i.c.w. CODE 24 side panel*	592,00
❑	781	USB charging socket	197,00
❑	782	Universal adapter for mounting mobile phone holders ⚠ *not in conjunction with R-NET OMNI II CODE 702/839 Cell phone holder not included in scope of delivery*	116,00
❑	4795	Anti-theft function, key ⚠ *only i.c.w. R-NET CODE 702*	54,00
❑	4772	TUV-certificate only i.c.w. code 117	0,00
❑	2721	Adapter set for DAHL docking station ⚠ **Rollstuhlseitiges Adapterset für die Vorbereitung der Dahl-Dockingstation* *Nicht i.V. mit Lenkradfederung für 8" Rad CODE 2145 und 2146*	1.147,00

Wheelchair options

	CODE	DESCRIPTION	EUR €
❑	2676	Warning sign sticker ⚠ *only in combination with ErgoSeat*	109,00
❑	676	Rear marking plate	109,00

EPOWER

☐ Please select by ticking the relevant boxes
Standard version: If no other selection is made, the standard version will be delivered.
S Special construction/individual solution o.r. - on request

24 Service

	CODE	Item	DESCRIPTION	EUR €
❑		55704	Flat fee for customised parameterisation (not eligible for discounts) ⚠ *The on-the-spot service covers only simple adjustment of parameter sets, such as acceleration, speed etc.*	556,00
❑		55709	Medically required parameterisation of the standing and lie-down module to customer-specific settings (not eligible for discounts) ⚠ *The on-the-spot service covers simple program adjustments as well as the integration of special controls.*	1.124,00