EUROCHAIR

1.750 Stock version

Allowed user weight

130 kg

Suitable for passenger transport

Tested to ISO 7176-19



MEYRA GROUP



Angle adjustable steering head
Aluminium steering fork



High variability through retrofit kits



Conversion possible quickly and specifically without exchanging parts



Legrests can be swivelled inwards and outwards, choice of 2 frame colours

For universal application

EFFICIENCY

- Durable product through the use of high value materials and a high level of manufacturing quality, little need for repairs
- Compatibility with MEYRA's modular system
- Long term assurance of spare parts provision through an area wide dealer structure
- Optimal reinstallment through simple adaptation possibilities to the user
- All basic settings are possible without exchanging parts

THERAPEUTIC BENEFIT

- Enhances independent mobility through fine tuned adaptation to the user
- High variability through retrofit kits
- Better light-running properties with little required propelling strength
- Simple handling and reduced weight ease everyday handling

FIELD OF USE

- For permanent everyday use
- For optimal provision of many disabilities due to the high level of adjustment possibilities
- Mostly implemented at home, but also in rehabilitation clinics to achieve therapeutic goals
- For big and tall people

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

MEYRA



Standard configuration

- Rear wheels 24"
- Handrim aluminium
- Handrim distance wide (can be changed to narrow)
- PU-tyres, puncture safe
- Steering wheel 150 mm PU
- Backrest with lumbar support
- Side guard can be folded up, swivelled away and removed
- Armrest pad depth adjustable
- Leg supports, removable and swivelling
- Seat depth 430 mm (adjustable)
- Seat and back strap Polyester, black
- Frame colour Silverline or Navy blue
- Adjustable back with back pocket
- Angle adjustable footplates
- Seat height front 500 mm
- Seat height rear 470 mm (adjustable)
- Back height 420 mm (adjustable)

Allowed user weight

130 kg

Suitable for passenger transport

Tested to ISO 7176-19

Hints

All prices are subject to the statutory value added tax!
Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
You can find information on the measurement process at www.meyra.de/messverfahren

Valid from 01.01.2026 - 31.12.2026.
Errors, misprints and changes reserved.

When transported in public transport (e.g. train), the dimensions of the wheelchair may be higher than the guide values of the transport company. Please inquire in advance.
When using auxiliary drives, the permissible total weight must not be exceeded.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

MANUAL

Technical specifications

[1] Applies to wheelchair with pressure brake

[2] Applies to wheelchairs with narrow handrim spacing and wheelbase

Seat widths in mm	380 / 400 / 430 / 460 / 480 / 500 / 530
Ready for use width in mm	560 / 580 / 610 / 640 / 660 / 680 / 710 [1,2]
Width folded in mm	280
Seat depth in mm	400 - 460
Seat height front/rear in mm	420 - 530 / 430 - 530
Overall height in mm	960
Back support height in mm	400 / 420
Arm support height in mm	from 230
Length with footrest in mm	from 1000
User weight	130 kg
Empty weight from	15 kg [1]
Permitted total weight in kg	150 kg
Rear wheel	600 (24")
Brake test according to STANDARD	7°

Please enter your quantities in the boxes

01 Pressure brake, frame colour silverline

Item	SW	DESCRIPTION	Ø
9173893	380 mm	Pressure brake, frame colour silverline	24"
9173894	400 mm	Pressure brake, frame colour silverline	24"
9173895	430 mm	Pressure brake, frame colour silverline	24"
9173896	460 mm	Pressure brake, frame colour silverline	24"
9173897	480 mm	Pressure brake, frame colour silverline	24"
9173898	500 mm	Pressure brake, frame colour silverline	24"
9173899	530 mm	Pressure brake, frame colour silverline	24"

02 Drum brake, frame colour silverline

Item	SW	DESCRIPTION	Ø
9173900	380 mm	Drum brake, frame colour silverline	24"
9173901	400 mm	Drum brake, frame colour silverline	24"
9173902	430 mm	Drum brake, frame colour silverline	24"
9173903	460 mm	Drum brake, frame colour silverline	24"
9173904	480 mm	Drum brake, frame colour silverline	24"
9173905	500 mm	Drum brake, frame colour silverline	24"
9173906	530 mm	Drum brake, frame colour silverline	24"

03 Pressure brake, frame colour navy blue

The wheelchair in the frame colour navy blue can only be ordered in the basic version as a direct container (270 units)!

Item	SW	DESCRIPTION	Ø
9173907	380 mm	Pressure brake, frame colour navy blue	24"
9173908	400 mm	Pressure brake, frame colour navy blue	24"
9173909	430 mm	Pressure brake, frame colour navy blue	24"
9173910	460 mm	Pressure brake, frame colour navy blue	24"
9173911	480 mm	Pressure brake, frame colour navy blue	24"
9173912	500 mm	Pressure brake, frame colour navy blue	24"
9173913	530 mm	Pressure brake, frame colour navy blue	24"

Please enter your quantities in the boxes

Retrofit sets

The following retrofit kits are enclosed and must be fitted by a specialist dealer!
For more information on fitting, see the service instructions!

04 Backrests

Item	DESCRIPTION
1090825	Back tube without lumbar bend, BH 400 mm / 420 mm
1090826	Back tube without lumbar bend with telescopic push handle 100 mm height adjustable, BH 420 mm
1090827	Back tube with lumbar bend with telescopic push handle 100 mm height adjustable, BH 420 mm
1090828	Back tube 10° angle adjustable with telescopic push handle 100 mm height adjustable, BH 420 mm

05 Arm supports

Item	DESCRIPTION
1090829	Side guard incl. height adjustable and removable armrest via push button, pair, RSH 470 mm

06 Leg supports

Item	DESCRIPTION
1090830	Leg support upper part, removable, swivel away and up, pair

07 Seat heights

Item	DESCRIPTION
1098194	Front seat height 420 - 500 mm, steering wheel 125 mm solid rubber
1098192	Front seat height 430 - 510 mm, steering wheel 142 mm soft solid rubber
1098193	Front seat height 485 - 530 mm, steering wheel 180 mm pneumatic tyres
1090832	Front seat height 485 - 530 mm, steering wheel 180 mm PU

Please enter your quantities in the boxes

08 Seat Cushions

Item	SW-SD-CH (mm)	DESCRIPTION
1046377	380-400-30	Seat cushion black
1006685	380-430-30	Seat cushion black
1006705	380-460-30	Seat cushion black
1046377	400-400-30	Seat cushion black
1006685	400-430-30	Seat cushion black
1006705	400-460-30	Seat cushion black
1006685	430-400-30	Seat cushion black
1006825	430-430-30	Seat cushion black
1015928	430-460-30	Seat cushion black
1006705	460-400-30	Seat cushion black
1015928	460-430-30	Seat cushion black
1015929	460-460-30	Seat cushion black
1046378	480-400-30	Seat cushion black
1046382	480-430-30	Seat cushion black
1046381	480-460-30	Seat cushion black
1046378	500-400-30	Seat cushion black
1046382	500-430-30	Seat cushion black
1046381	500-460-30	Seat cushion black
1049053	530-400-30	Seat cushion black
1049054	530-430-30	Seat cushion black
1049055	530-460-30	Seat cushion black

Please enter your quantities in the boxes

09 Other accessories

Item	DESCRIPTION
1041414	Safety belt (metal buckle)
1063168	Crutch holder
1047486	Reflectors
1090833	Brake lever extension attachable, pair
1090835	Anti tip wheels, adjustable length, removable, pair
1069874	Safety package, consisting of passive lighting, support castors, retaining strap with lock, walking aid holder