FORMAT

3.940 Stock version

Allowed user weight

130 kg

Suitable for passenger transport

Tested to ISO 7176-19



ORTOPEDIA

EINE MARKE VON MEYRA




Height-adjustable side panels as standard



Height-adjustable push handles as standard



High-quality adjustable backrest as standard



Angle adjustable footplates

The quality alternative for reinstallment

EFFICIENCY

- The economic start into the class of lightweight wheelchairs
- Due to the dependence from the proven model Service 3.600 and Eurochair 1.750 the 3.940 guarantees perfect and uncomplicated spare parts provision from the already existing product program

THERAPEUTIC BENEFIT

- In case of increased adaptation demands, the Format supports the user when standing up and transferring, since the seat height can be aligned to the body size of the user.

FIELD OF USE

- Format is perfectly suited for a start into your personal mobility
- Push- handrim drive
- Home and inpatient application
- Rheumatism
- Geriatric provision

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

ORTOPEDIA
EINE MARKE VON MEYRA



Standard configuration

- Drive wheel 24" aluminium
- Steering wheel 180 x 50 mm
- Puncture-proof tyres
- Adjustable back with back pocket
- Backrest with lumbar support
- Push handles 100 mm height-adjustable, telescopic
- Side part height-adjustable, fold-up-, swivable-, detachable
- Leg supports, removable and swivelling
- Foot plates angle adjustable
- Seat heiht 510 mm
- Seat depth 430 mm
- Back height 420 mm (adjustable)
- Handrim distance narrow
- Frame colour Silverline

Allowed user weight

130 kg

Suitable for passenger transport

Tested to ISO 7176-19

Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
When transported in public transport (e.g. train), the dimensions of the wheelchair may be higher than the guide values of the transport company. Please inquire in advance.
When using auxiliary drives, the permissible total weight must not be exceeded.

Valid from 01.01.2026 - 31.12.2026.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

MANUAL

Technical specifications

Seat widths in mm	380 / 400 / 430 / 460 / 480 / 510
Width ready to drive in mm from	Seat width + 180
Width folded in mm from	280
Seat depth in mm	400 / 430
Front seat height in mm	470 - 510
Seat height in the back in mm	440 - 490
Overall height in mm	960
Back support height in mm	420
Arm support height in mm	215
Length with footrest in mm	1.040
User weight	130 kg
Empty weight from	17,5 kg
Transport weight	10,5 kg
Rear wheel	560 mm (22") / 610 mm (24")
Brake test according to STANDARD	7°

Please enter your quantities in the boxes

01 Pressure brake (push to lock)

Item	DESCRIPTION
9394025	SW 380 mm Pressure brake
9394026	SW 400 mm Pressure Brake
9394027	SW 430 mm Pressure Brake
9394028	SW 460 mm Pressure brake
9394029	SW 480 mm Pressure brake
9394030	SW 510 mm Pressure brake

02 Drum Brake

Item	DESCRIPTION
9394045	SW 380 mm Drum brake
9394046	SW 400 mm Drum brake
9394047	SW 430 mm Drum brake
9394048	SW 460 mm Drum brake
9394049	SW 480 mm Drum brake
9394050	SW 510 mm Drum brake

03 Options (delivered separately)

Item	DESCRIPTION
1069874	Safety package, consisting of: Retaining strap, cane holder, passive lighting, support castors
1064946	Therapy table plastic, transparent, adjustable from 380 - 480 mm seat width

Please enter your quantities in the boxes

04 Seat cushion (supplied separately)

Item	SW-SD-CH (mm)	DESCRIPTION
1046377	380-400-30	Seat cushion black
1006685	380-430-30	Seat cushion black
1046377	400-400-30	Seat cushion black
1006685	400-430-30	Seat cushion black
1006685	430-400-30	Seat cushion black
1006825	430-430-30	Seat cushion black
1006705	450-400-30	Seat cushion black
1006825	450-430-30	Seat cushion black
1046378	480-400-30	Seat cushion black
1046382	480-430-30	Seat cushion black