HM303D HM03MS

7.010 Stock version



MEYRA GROUP



HM303D



HM03MS

SUPER LIGHT WHEELCHAIR

CHARACTERISTICS

- This wheelchair is designed for older people to gain more flexibility in everyday life, for more independence
- A compact wheelchair that is very easy to carry and load
- For more safety in everyday life
- This wheelchair is so light that it is very easy to lift into the car, it folds easily and is very manoeuvrable
- Drive wheels with quick-release axle for space-saving storage
- Swivelling leg rests for more freedom of movement

● ORDER ❍ COST ESTIMATE ❍ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

MEYRA



Standard configuration

- HM303D
- Frame colour matt black
- Low weight - only 10 kg
- Foldable back
- pelvic belt
- Dual brake system
- Transport wheelchair
- Portable
- Drive wheels 22'

- HM03MS
- Frame colour navy blue
- Low weight - only 9,6 kg
- Foldable back
- pelvic belt
- Dual brake system
- Transport wheelchair
- Portable
- Drive wheels 12'

Hints

All prices are subject to the statutory value added tax!
Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
Valid from 01.01.2026 – 31.12.2026.
Errors, misprints and changes reserved.
Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

MANUAL

Technical specifications

	HM303D	HM03MS
Seat width in mm	450	450
Seat depth in mm	410	410
Seat height in mm	450	465
Back support height in mm	400	410
Backrest Angle	82° +/- 1°	82° +/- 1°
Arm support height in mm	220	220
Ready for use width in mm	645	600
Overall height in mm	860	875
Total weight in kg	10 kg	9,6 kg
Footrest length in mm	390	390
Overall length with footrest in mm	942	1075
Max. User weight in kg	110 kg	110 kg
PU wheels	6"/22"	6"/12"
Packaging dimensions in mm	930 x 360 x 710	910 x 350 x 720

Version No.: #1554 - 16.04.2026 / 2026-08-25 11:04:21

Please enter your quantities in the boxes

HM303D

Item	DESCRIPTION	EUR €
9350822	HM303D 22" SB 450 - ST 410 - SH 450 DUAL BRAKE SYSTEM	814,00

HM03MS

Item	DESCRIPTION	EUR €
9350812	HM03MS 12" SB 450 - ST 410 - SH 465 DUAL BRAKE SYSTEM	710,00