BUDGET

9.050 Stock version

Allowed user weight

130 kg

Suitable for passenger transport

Tested to ISO 7176-19



ORTOPEDIA

EINE MARKE VON MEYRA



Easy removal function of the side plate



Knee lever brake



Lockable armrests can be folded back



Front seat height adjustable without changing parts

Economic and innovative

EFFICIENCY

- Specially designed for the current provision situation
- Technical variability through compatibility with EUROCHAIR: Arm supports, leg supports, seat- and back straps, pressure and drum brake

THERAPEUTIC BENEFIT

- Easy to fold for easy handling in everyday life

FIELD OF USE

- Push- handrim drive
- Home and inpatient application
- Rheumatism
- Geriatric provision

● ORDER ○ COST ESTIMATE ○ CUSTOM DESIGN

Please send the order form by e-mail to: info@meyragroup.com **Customer service:** +49 571 93292 555

Billing address / Customer ref. no.:

Delivery address:

Tax-ID

Consignment / Comments:

ORTOPEDIA
EINE MARKE VON MEYRA



Standard configuration

- Drive wheel 24" aluminium
- Steering wheel 180 x 50 mm
- Puncture-proof tyres
- Arm supports up-foldable
- Leg supports, removable and swivelling
- Seat heiht 510 mm
- Seat depth 430 mm
- Back height 420 mm
- Handrim distance narrow
- Frame colour Silverline
- Backrest with lumbar support

Allowed user weight

130 kg

Suitable for passenger transport

Tested to ISO 7176-19

Hints

Our Terms and Conditions as stated at www.meyra.de/AGB shall apply
When transported in public transport (e.g. train), the dimensions of the wheelchair may be higher than the guide values of the transport company. Please inquire in advance.
When using auxiliary drives, the permissible total weight must not be exceeded.

Valid from 01.01.2026 – 31.12.2026.
Errors, misprints and changes reserved.

Use our K.I.S.S.S. for orders and configurations at www.kisss-by-meyra.de
You avoid queries and incorrect orders thanks to the buildability check

MANUAL

Technical specifications

Seat widths in mm	380 / 400 / 430 / 460 / 480 / 510
Width ready to drive in mm from	Seat width + 180
Width folded in mm from	280
Seat depth in mm	430
Front seat height in mm	470 - 510
Seat height in the back in mm	440 - 490
Overall height in mm	960
Back support height in mm	420
Arm support height in mm	210,5
Length with footrest in mm	1040
User weight	130 kg
Empty weight from	19,6 kg
Permitted total weight in kg	150 kg
Rear wheel	610 mm (24")
Brake test according to STANDARD	7°

MANUAL

Please enter your quantities in the boxes

01 Pressure brake (push to lock)

Item	DESCRIPTION
9905141	SW 380 mm Pressure brake
9905142	SW 400 mm Pressure Brake
9905143	SW 430 mm Pressure Brake
9905144	SW 460 mm Pressure brake
9905145	SW 480 mm Pressure brake
9905146	SW 510 mm Pressure brake

02 Drum Brake

Item	DESCRIPTION
9905147	SW 380 mm Drum brake
9905148	SW 400 mm Drum brake
9905149	SW 430 mm Drum brake
9905150	SW 460 mm Drum brake
9905151	SW 480 mm Drum brake
9905152	SW 510 mm Drum brake

Please enter your quantities in the boxes

03 Options (delivered separately)

Item	DESCRIPTION
1069874	Safety package, consisting of: Retaining strap, cane holder, passive lighting, support castors

04 Seat cushion (supplied separately)

Item	SW-SD-CH (mm)	DESCRIPTION
1006685	380-430-30	Seat cushion black
1006685	400-430-30	Seat cushion black
1006825	430-430-30	Seat cushion black
1015928	460-430-30	Seat cushion black
1046382	480-430-30	Seat cushion black
1046382	510-430-30	Seat cushion black